





CLARITY® LC SOLUTIONS

Partner in Purification

The number of requests for high purity oligonucleotides with synthesis scales in the µmole range (and greater) is growing and has created a need for more efficient and higher capacity purification solutions. Phenomenex would like to be your partner in purification and assist you in supplying small to large quantities of high purity synthetic DNA/RNA to your valued customers. As your purification partner, we not only provide excellent technical support and customer service, but also a novel, reliable suite of products (Clarity BioSolutions) for synthetic DNA/RNA purification. We are pleased to offer a line of LC purification solutions within the Clarity BioSolutions portfolio. Clarity Oligo-RP™ LC columns are an excellent product for those who require a high purity reversed phase solution with long lifetime, while Clarity Oligo-WAX™ LC columns are suitable for those who require an ion exchange solution with high capacity. We look forward to being a partner to companies and core labs who demand efficient, economical, and efficacious synthetic oligo purification and support.

Table of Contents

| F Clarity® Oligo-RP™: | PAGE |
|--|-------|
| High Purity, High Efficiency Reversed Phase LC Solution | 1 - 6 |
| Clarity® Oligo-WAX™: High Purity, High Loadability Ion Exchange LC Solution | 7 - 9 |
| Material Characteristics | . 10 |
| Ordering Information | . 10 |
| Clarity BioSolutions Portfolio | 11 |




High Purity, High Efficiency Reversed Phase LC Solution

larity Oligo-RP LC columns have been specifically designed for the reversed phase purification of oligonucleotides with balanced hydrophobicity and polar selectivity. The media is based on composite particle TWIN[™] technology. This technology gives improved selectivity and efficiency for oligonucleotides when compared to other hybrid, polymer, and silica particles found in the marketplace. It is available in 3, 5, and 10 µm particle sized beads, in a variety of dimensions, and utilizes HPLC instrumentation for increased productivity.

- Easily separate N-1 failure sequences from target oligo with > 90 % purities
- Trityl-off purification of DNA, RNA, thioates, and modified/labeled oligonucleotides
- Preparative dimensions & particle sizes for loads > 5 µmole
- Purify oligos up to 60 nt in length
- **Excellent column for reversed phase HPLC quality control (QC) testing**

RNA Purification of Failure N-1 Sequence from Target N Sequence



Column: Clarity 3 um Oligo-RP C18

Dimensions: 50 x 4 6 mm Part No.: 00B-4441-E0 Mobile Phase: A: 10 mM TPAC B: Methanol

Gradient: A/B (75:25) to A/B (55:45) in 80 min (50 min run time)


Flow Rate: 1 mL/min UV @ 260 nm Detection:

1. 20 nt RNA with sequence Sample:

CUGUAAUCUCUUGUCUATT (2.5 µg) 2. 21 nt RNA with sequence UCUGUAAUCUCUUGUCUATT (2.5 μg)

Excellent selectivity characteristics of Clarity Oligo-RP allow baseline separation of failure N-1 sequences from target N sequences.

Resolution Achieved Between 60 nt Impurities with Similar Sequences



Column: Clarity 3 um Oligo-RP C18

Dimensions: 50 x 4.6 mm 00B-4441-E0 Part No.:

Mobile Phase: A: 50 mM TEAA pH 7.5 / 5 % Acetonitrile

Gradient: 20 % to 25 % B in 20 minutes; hold at 5 minutes

@ 25 % B

Flow Rate: 1 mL/min UV @ 260 nm **Detection:**

60 nt DNA with sequence 5'-CTC CTG GGC CGT Sample:

GGC TCT GCG CAC TTC AGG AAA CTG GGC ACT

CCT GGG CAG TGG ATC TGC-3'

The high efficiency and selectivity of the sorbent as well as ion-pairing interactions produce a fingerprint of a crude 60 nt DNA on Clarity Oligo-RP illustrating resolution of impurities in the final product. Oligo-RP can recognize even the slightest changes in a nucleotide sequence. (numbers on each peak represent the sequence length of the impurity present)

High Purity, High Efficiency Reversed Phase LC Solution



TWIN[™] Technology – Engineered for DNA/RNA Purification


The Oligo-RP media is based on the TWIN technology used to manufacture the patent pending, state-of-the-art Gemini® material, but has been tailored to purify synthetic, as well as natural, DNA & RNA based on the chemical characteristics of these molecules. This solution is able to recognize minute differences in interaction features of two oligonucleotides with very closely resembling structures (for example N/N-1 sequences). In addition, this recognition ability enables Clarity Oligo-RP to provide better resolution between such close pairs of sequences both on the analytical and preparative scale.

- Long column lifetime due to extended pH stability & mechanical strength
- Excellent efficiency enables separation of oligos with similar chemistries
- Improved oligo selectivity over hybrid, polymer, and silica particles

Manufacturing of TWIN Media

TWIN[™] **Technology**


During the final stage of silica manufacturing a unique silica-organic layer is grafted to create a completely new composite particle. Since the internal base silica is unaltered by this manufacturing process, the particle retains its mechanical strength and rigidity along with excellent efficiency, while the silica-organic shell protects the particle from chemical attack.



High Purity, High Efficiency Reversed Phase LC Solution

TWIN™ Technology – Engineered for DNA/RNA Purification (cont'd)

Extended Lifetime and High Efficiency



Mobile Phase: A: 10 mM Ammonium Bicarbonate, pH 10.0


B: 90:10 Acetonitrile/Buffer

Gradient: 0 % to 100 % B in 10 min; hold at 100 % B for 7 min; re-equilibrate at 0 % B for 3 min

1.0 mL/min Flow Rate: Temperature: 50 °C Detection: UV @ 254 nm 1. Amitryptyline Sample: 2. Prednisolone

The unique engineering of Clarity Oligo-RP provides stability, high efficiency, and increased column lifetime compared to other commercial **HPLC** columns.

TWIN vs Hybrid Technology



High Purity, High Efficiency Reversed Phase LC Solution




Advantages of Utilizing Reversed Phase HPLC


eversed phase separation of oligonucleotides has advantages over other modes of separations such as ion exchange and PAGE. The Oligo-RP phase allows high loadability and delivers high recovery and purity, eliminating the need for extra purification steps such as desalting. Because reversed phase HPLC is an extremely high-resolution technique, Oligo-RP columns are easily able to separate the target oligonucleotide from the many impurities in the mixture, such as the N-1 failure sequence. In addition, this one phase can separate both DNA and RNA chemistries and modified oligos.

- Recognize minute differences in oligonucleotides such as N and N-1 sequences
- Desalting not required after HPLC analysis
- Base line separation/collection of desired peak at preparative scale (5 μmole and greater)
- Compatible with MS friendly buffers
- 3 µm material available for improved efficiency and resolution

High Purity Preparative DNA Purification



40 nt DNA Purification from 39 nt Impurity




The combination of reversed phase conditions with an Oligo-RP column can recognize even the slightest changes in a nucleotide sequence, such as a substitution of one base by another or a difference of one base as seen in this application.

High Purity, High Efficiency Reversed Phase LC Solution

Clarity Oligo-RP Reversed Phase Applications

ligo-RP can be used for a variety of oligonucleotide analyses and purifications including RNA, DNA, modified oligos, labeled sequences, and many more. If you have specific questions on whether Oligo-RP is the right product for you, please contact Phenomenex and we will assist you in selecting the correct solution.

Separate ssDNA from dsDNA



Column: Clarity 3 um Oligo-RP **Dimensions:** 50 x 4.6 mm ID Part No.: 00B-4441-E0 Mobile Phase: A: 50 mM TEAA B: MeOH

Gradient: A/B (90:10) to A/B (40:60) in 20 min Flow Rate: 1 mL/min


Detection:

UV-Vis Abs. - Variable Wave.

Temperature: ambient Sample: 1. ss DNA 2. ds DNA

Clarity Oligo-RP is capable of baseline separating single stranded from double stranded DNA due to its unique selectivity.

RP-HPLC Purification of 40 nt Oligos with Varying GC Content



Column: Clarity 3 µm Oligo-RF **Dimensions:** 50 x 4.6 mm ID Part No.: 00B-4441-E0

Mobile Phase: A: 50 mM TEAA / 5 % Acetonitrile

B: Methanol

Gradient: A/B (90:10) to A/B (40:60) in 20 min

Flow Rate: 1 mL/min UV @ 260 nm Detection: Temperature: ambient 1. DNA 40 nt (GC: 65 %)

Sample:


2. DNA 40 nt (GC: 50 %) 3. DNA 40 nt (GC: 37 %)

C & G rich compounds were separated from non C & G rich ones. This was possible due to the fact that the sorbent possesses a tailored mixture of hydrophobic, dipolar, pi-pi, and hydrogen bonding characteristics.

High Purity, High Efficiency Reversed Phase LC Solution



Dye-labeled DNA Purification



 Column:
 Clarity 3 µm Oligo-RP

 Dimensions:
 50 x 4.6 mm ID

 Part No.:
 00B-4441-E0

 Mobile Phase:
 A: 50 mM TEAA / 5 % Acetonitrile

 B: Methanol

 Gradient:
 A/B (90:10) to A/B (60:40) in 20 min

Flow Rate: 1 mL/min
Detection: UV @ 260 nm


Sample: 1. FAM-DNA (25 nt) GC% 72 %

ambient

2. FAM Labeled DNA (25 nt) GC% 56 % 3. FAM Labeled DNA (25 nt) GC% 20 %

Oligo-RP delivers high purity labeled oligos such as those with a FAM label. Baseline separation is achieved ensuring good yields and purity.

40 nt DNA Purification



 Column:
 Clarity 3 µm Oligo-RP

 Dimensions:
 50 x 4.6 mm ID

 Part No.:
 00B-4441-E0


Mobile Phase: A: 50 mM TEAA / 5 % Acetonitrile B: Methanol

Gradient: A/B (90:10) to A/B (40:60) in 20 min

Flow Rate: 1 mL/min
Detection: UV @ 260 nm
Temperature: ambient
Sample: 1. DNA 40 nt

40 nt purifications are sometimes difficult on a reversed phase column. Oligo-RP, however, can purify lengths up to 60 nt and sometimes longer depending on the sequence.

21 nt RNA Purification



 Column:
 Clarity 3 μm Oligo-RP

 Dimensions:
 50 x 4.6 mm ID

 Part No.:
 00B-4441-F0

Mobile Phase: A: 50 mM TEAA / 5 % Acetonitrile

B: Methanol

Gradient: A/B (90:10) to A/B (40:60) in 20 min

 Flow Rate:
 1 mL/min

 Detection:
 UV @ 260 nm

 Temperature:
 ambient

 Sample:
 1. RNA 21 nt

RNA, due to the fact that it is more polar, is more of a challenge to separate on a RP-HPLC columns than DNA. The dipolar and hydrogen binding aspect of Oligo-RP makes it an excellent purification tool for RNA sequences such as the 21 nt in this application.


CLARITY® OLIGO-WAX™

High Purity, High Loadability Preparative Ion Exchange LC Solution

larity Oligo-WAX LC columns were designed with the synthetic DNA/RNA preparative chromatographer in mind. In preparative chromatography, it is imperative that not only high purities are achieved but also that the media has high capacity and efficiency to produce high yields. Oligo-WAX is an advantageous combination of purity, capacity, mechanical strength, and efficiency.

- Columns amenable to HPLC and FPLC systems
- High capacity media
- Increased efficiency versus other IEX media
- High purity silica-based media

20 nt DNA Purification



Column:Clarity 10 μm Oligo-WAXDimensions:150 x 4.6 mm ID

Part No.: 00F-4451-E0

Mobile Phase: A: 20 mM Tris pH 8.0 / 10 % Acetonitrile B: 20 mM Tris pH 8.0 / 10 % Acetonitrile /

1.2 M NaCl

Gradient: 100 % A to 100 % B in 20 min

Flow Rate: 2.2 mL/min
Detection: UV-Vis Abs.

Detection: UV-Vis Abs. - Diode Array 260 nm **Temperature:** ambient

Sample: 1. Desalted DNA 20 nt

After desalting, a 75 µg mixture of a synthesized 20 nt ssDNA (ACG TCA TGT CGA GAT CAT CG) was loaded onto a Clarity Oligo-WAX column to purify it from closely eluting impurities. Note the good separation, due to the high efficiency and selectivity of the silica-based weak anion exchanger, of the full-length oligo from N-1 impurities and other contaminants.

CLARITY® OLIGO-WAX™






Oligo-WAX Technology – Tailored for Preparative Purification

he majority of synthetic oligo preparative purifications are performed using a strong anion exchanger bonded to a $10 \text{ or } 15 \mu\text{m}$ polymer backbone. Polymer backbones are amenable to clean in place protocols, while strong anion exchangers have a wide effective pH range. To date, these technologies have been satisfactory for prep purifications and will continue to be. However, due to the fact that Clarity Oligo-WAX is a crosslinked weak anion exchanger bonded to a $10 \mu\text{m}$ high purity silica, this technology offers a few advantages that aren't available with typically used purification products.

- Higher loading capacity due to very high density ligand
- Improved peak efficiency & flow characteristics due to pore morphology and spherical shape of the silica matrix
- Better fractionation of closely eluting compounds resulting in higher purity
- Increased productivity by running at higher flower rates and pressures

CE Purity Analysis of Ion Exchange Purification



- Final Purity = 88.8 %, N-1 = 2.0 % Final amount = 205.9 OD's
- Recovery of full-length product = 28.9 % Conductivity = 200 μ S/cm
- Final Purity = 95.1 %, N-1 = 1.3 % Final amount = 188.9 OD's

• Recovery of full-length product = 28.4 % • Conductivity = 151 μS/cm

Two purification runs were performed on each column with fractional QC being taken after each run. Passing fractions from the two purification runs were combined into one pooled lot for each column. That pooled lot was then divided equally and run through a Clarity desalting tube. Final OD's and QC were taken after desalting, including ESI, CE, and conductivity. The purity and resolution of Clarity Oligo-WAX was considerably better than SOURCE 15Q. Though SOURCE had a slightly higher recovery of full length oligo, it was not be a wide enough margin to offset the purity advantage.

Data courtesy of a large, lowa-based oligo manufacturer. Comparative separations may not be representative of all applications.


CLARITY® OLIGO-WAX"

High Purity, High Loadability Preparative Ion Exchange LC Solution

Clarity Oligo-WAX Ion Exchange Applications

ligo-WAX is recommended for preparative purifications of synthetic RNA and DNA. If you have specific questions on whether Oligo-WAX is the right product for you or how it differs from commercial SAX columns on the market, please contact Phenomenex and we will be pleased to assist you.

Purification of Long Oligo Sequences



 Column:
 Clarity 10 μm Oligo-WAX

 Dimensions:
 150 x 4.6 mm ID

 Part No.:
 00F-4451-E0

Mobile Phase: A: 20 mM Tris pH 8.0 / 10 % Acetonitrile

B: 20 mM Tris pH 8.0 / 10 % Acetonitrile / 1.2 M NaCl

100 % A to 100 % B in 20 min

Flow Rate: 2.2 mL/min

Detection: UV-Vis Abs. - Diode Array 260 nm


Temperature: ambient

Gradient:

Sample: 1. DNA Purified 40 nt

100 µg for a 40 nt ssDNA (CTC CTG GGC AGT GGA TCT GCG CAC TTC AGG CTC CTG GCG A) was purified using Clarity Oligo-WAX. In purifying longer sequences, such as this 40 nt DNA, the results produced from commonly used ion exchange columns are typically broad peaks. However, the efficiency of the Oligo-WAX column generates sharp peaks, which are excellent for better purification of closely eluting contaminants.

DNA Purification of Hydrolyzed N-1 from N Sequence



 Column:
 Clarity 10 µm Oligo-WAX

 Dimensions:
 150 x 4.6 mm ID

 Part No.:
 00F-4451-E0

Mobile Phase: A: 20 mM Tris pH 8.0 / 10 % Acetonitrile B: 20 mM Tris pH 8.0 / 10 % Acetonitrile /

1.2 M NaCl

Gradient: 100 % A to 100 % B in 20 min

Flow Rate: 2.2 mL/min

Detection: UV-Vis Abs. - Diode Array 260 nm **Sample:** 1. Depurinated A & G & 20 nt DNA

To test the separation efficiency of the Oligo-WAX column, a mixture of a 20 nt synthesized DNA (ACTCGGCTTCCTCCTCTT) was exposed to polyamine hydrolysis to generate a 19 nt contaminant. This mixture was run on the Oligo-WAX column. Note the near baseline separation between the almost identical 19 nt and 20 nt.

CLARITY® LC SOLUTIONS

Reversed phase & ion exchange high purity synthetic oligonucleotide purification

Material Characteristics

| | Clarity® Oligo-RP™ Clarity® Oligo-WAX™ | |
|--------------------|--|---|
| Particle Support | TWIN [™] (silica organic composite) | High Purity Silica |
| Bonded Phase | C18 | Crosslinked polyamine (WAX) |
| Particle Size | 3, 5, and 10 µm | 10 μm |
| Pore Size | 110 Å | 360 Å |
| Surface Area | 375 m²/g | 160 m²/g |
| pH stability | 1.0 – 12.0 | 1.5 – 11.0 |
| pH Effective Range | 1.0 – 12.0 | 4.0 - 8.5 (Amine ligand deprotonates at pH 9) |

Ordering Information

| 3 µm Columns (mm) | | | | | |
|-------------------|-------------|-------------|-------------|-------------|-------------|
| Phases | 50 x 1.0 | 100 x 1.0 | 50 x 2.0 | 100 x 2.0 | 150 x 2.0 |
| Oligo-RP | 00B-4441-A0 | 00D-4441-A0 | 00B-4441-B0 | 00D-4441-B0 | 00F-4441-B0 |
| Phases | 50 x 4.6 | 100 x 4.6 | 150 x 4.6 | 50 x 10.0 | |
| Oligo-RP | 00B-4441-E0 | 00D-4441-E0 | 00F-4441-E0 | 00B-4441-N0 | |

| 5 μm Columns (mm) | | | | | |
|-------------------|-------------|-------------|-------------|-------------|-------------|
| Phases | 50 x 1.0 | 150 x 2.0 | 50 x 4.6 | 150 x 4.6 | 250 x 4.6 |
| Oligo-RP | 00B-4442-A0 | 00F-4442-B0 | 00B-4442-E0 | 00F-4442-E0 | 00G-4442-E0 |
| Phases | 50 x 10.0 | 100 x 10.0 | 250 x 10.0 | 250 x 21.2 | 250 x 30.0 |
| Oligo-RP | 00B-4442-N0 | 00D-4442-N0 | 00G-4442-N0 | 00G-4442-P0 | 00G-4442-U0 |

| 10 µm Columns (mm) | | | | |
|--------------------|-------------|-------------|-------------|-------------|
| Phases | 100 x 4.6 | 150 x 4.6 | 150 x 10.0 | 250 x 10.0 |
| Oligo-RP | _ | 00F-4445-E0 | 00F-4445-N0 | - |
| Oligo-WAX | 00D-4451-E0 | 00F-4451-E0 | 00F-4451-N0 | 00G-4451-N0 |
| Phases | 250 x 21.2 | 150 x 30.0 | 150 x 50.0 | 250 x 50.0 |
| Oligo-RP | 00G-4445-P0 | 00F-4445-U0 | 00F-4445-V0 | 00G-4445-V0 |
| Oligo-WAX | _ | 00F-4451-U0 | _ | - |

GUARD CARTRIDGE SYSTEM

| Part No. | Descri | Unit | For ID | |
|----------|---------------------------------------|---|--------|--------------|
| AJ0-8134 | SecurityGuard [™] Cartridges | Clarity Oligo-RP 4 x 2.0 mm* | 10/Pk | 2.0 - 3.0 mm |
| AJ0-8135 | SecurityGuard Cartridges | Clarity Oligo-RP 4 x 3.0 mm* | 10/Pk | 3.2 - 8.0 mm |
| AJ0-8136 | SecurityGuard Semi-Prep Cartridges | Clarity Oligo-RP 10 x 10 mm [‡] | 3/Pk | 9 - 16 mm |
| AJ0-8210 | SecurityGuard PREP Cartridge | Clarity Oligo-RP 15 x 21.2 mm ID** | Ea | 18 - 29 mm |
| AJ0-8310 | SecurityGuard PREP Cartridge | Clarity Oligo-RP 15 x 30 mm ID** | Ea | 30 - 49 mm |
| AJ0-8324 | SecurityGuard Cartridges | Clarity Oligo-WAX 4 x 3.0 mm* | 10/Pk | 3.2 - 8.0 mm |
| AJ0-8325 | SecurityGuard Semi-Prep Cartridges | Clarity Oligo-WAX 10 x 10.0 mm [‡] | 3/Pk | 9 - 16 mm |

^{*}SecurityGuard Analytical Cartidges require holder, Part No.: KJ0-4282 ‡SecurityGuard Analytical Cartidges require holder, Part No.: AJ0-7220

^{**}SecurityGuard Analytical Cartidges require holder, Part No.: AJ0-8223

CLARITY LC SOLUTIONS

Clarity BioSolutions Portfolio

henomenex also offers other reliable products for synthetic DNA/RNA purification. Clarity® QSP™ is an excellent product for those who require a high-throughput, high purity solution. Clarity Desalting Tubes are a polyfunctional silica-based C18 sorbent that provides a high capacity, fast and effective desalting process. For more information, please visit www.phenomenex.com/clarity.

Clarity QSP 96-Well Plates

- For high-throughput, parallel processing of purified oligos
- Load 0.2 μmole synthesis scale per well
- Purify 96 crude oligo samples in ~45 minutes
- Easily amenable to automated liquid handling system

Clarity QSP Cartridge

- 50 mg/1 mL for purifying 0.2 μmole synthesis scale or less
- 150 mg/ 3 mL for purifying up to 1.0 µmole scale
- 5 g/ 60 mL for large scale purifications up to 50 μmole
- Purify crude oligo samples in ~8 minutes
- Use either vacuum or positive pressure systems

Clarity Desalting Tubes

- For the removal of salt and excess reagent
- Provides moderate purification of synthetic oligos









Evaluate Clarity® Biosolutions in your lab for 45 days, if you are not completely satisfied return it for a full refund.



PARTNER IN PURIFICATION SM



www.phenomenex.com

www.pnenomenex.com
Phenomenex products are available worldwide. For the distributor in your country,
contact Phenomenex USA, International Department by telephone, fax or email: international@phenomenex.com.

Australia 🐺 PO Box 4084

Lane Cove, NSW 2066 Australia

tel.: 02-9428-6444 fax: 02-9428-6445 email: info@

phenomenex.com.au

Ireland Queens Avenue, Hurdsfield Ind. Est. Macclesfield, Cheshire SK10 2BN, UK

tel.: 01 247 5405 fax: +44 1625-501796 email: eireinfo@ phenomenex.com 01-319-1301 01-319-1300

Austria

anfrage@ phenomenex.com

Zeppelinstr. 5 63741 Aschaffenburg Germany

Via Emilia, 51/C 40011 Anzola Emilia (BO) Italy

051 736176 051 735302 italiainfo@ phenomenex.com Canada 411 Madrid Ave.

Torrance, CA 90501-1430 USA

(800) 543-3681 (310) 328-7768 phenomenex.com

New Zealand P O Box 31-601 Milford 0741 North Shore City New Zealand

09-4780951 09-4780952 info@ phenomenex.co.nz Denmark Gydevang 39-41 3450 Allerød Denmark

> 4824 8048 4810 6265

dkinfo@ phenomenex.com

Puerto Rico 273 Sierra Morena, Suite #104 San Juan, Puerto Rico 00926

(800) 541-HPLC (310) 328-7768 info@ phenomenex.com France

Parc des Grillons, Bat.3 60 route de Sartrouville 78232 Le Pecq Cedex France

01 30 09 21 10 01 30 09 21 11 franceinfo@ phenomenex.com

United Kingdom

Queens Avenue, Hurdsfield Ind. Est. Macclesfield, Cheshire SK10 2BN, UK

01625-501367 01625-501796 ukinfo@ phenomenex.com **Germany** Zeppelinstr. 5

Ophenomenex

63741 Aschaffenburg Germany

.breaking with tradition

06021-58830-0 06021-58830-11 anfrage@ phenomenex.com

USA

411 Madrid Ave. Torrance, CA 90501-1430 USA

(310) 212-0555 (310) 328-7768 info@ phenomenex.com

CLARITY® LC SOLUTIONS
Reversed Phase & Ion Exchange Columns
for high purity synthetic oligonucleotide purification





