



Core-Shell Technology for Proteins and Peptides

Ultra-High Resolution and Performance on HPLC and UHPLC Systems



phenomenex[®]
...breaking with traditionSM



Welcome to the Future of BioSeparations

Introducing **Aeris™**, a specialized line of reversed phase core-shell HPLC / UHPLC columns, built exclusively for the ultra-high performance separation and analysis of proteins and peptides.

These columns can provide improved **resolving power, selectivity, throughput, sensitivity, column lifetime,** and **method flexibility** compared to other fully porous and core-shell columns typically used for bioseparations.

Choose your optimal
Aeris column

See page 6!

Aeris WIDEPORE

p16

Large pore optimized for
protein diffusion



XB-C18

XB-C8

C4

Multiple selectivities

3.6 μm particle for
HPLC and UHPLC systems

3.6 μm

Aeris PEPTIDE

p26

Small pore optimized
for peptide diffusion



XB-C18

Ideal surface chemistry
for resolving peptides

1.7 μm

Two scalable particle sizes
for method and system flexibility

3.6 μm



Table of Contents

Aeris Core-Shell Technology

- p. 4 Core-Shell Particles Precision Engineered for Protein and Peptide Separations
- p. 6 Easy Column Selection
- p. 8 Benefits of Using Aeris Core-Shell Columns
- p. 16 Aeris WIDEPORÉ for Intact Proteins and Polypeptides
- p. 26 Aeris PEPTIDE for Peptides and Peptide Mapping
- p. 34 Method Development and Optimization Services
- p. 35 Ordering Information

“Aeris had better separation abilities and lower backpressure than other core-shell reverse phase columns I've tried with proteins. I had very good technical and customer support throughout the entire process! I'm very glad to have switched to Phenomenex columns!”


-LYNN PRUISNER, TECHNOLOGY COMPANY

Core-Shell Particles Precision Engineered for Protein and Peptide Separations


Core-shell particle technology provides **striking increases in peak capacity and resolution** at lower backpressures, giving chromatographers the ability to achieve ultra-high performance on ANY system, HPLC or UHPLC.

A uniform porous silica layer is grown around a solid, spherical silica core, providing effective retention and selectivity with improved resolution, speed, and recovery. Next, optimizing the pore size and shell thickness for intact proteins or smaller peptide fragments provides well-defined depth penetration of biomolecules leading to **maximum separation power**.


TEM and SEM of Aeris™ PEPTIDE 3.6µm Core-Shell Particles



**Cross section of an
Aeris core-shell particle**



Magnified cross section of the porous “shell”




Uniform particle size and shape

The precise architecture of core-shell particles provides dramatic leaps in performance in two important ways:


1

The thin, porous layer, or “shell”, decreases the diffusion path length, thus reducing the time it takes for biomolecules to adsorb/desorb into and out of the particle.

Aeris Core-Shell Particle





Fully Porous Particle





2

Uniform sizing and shape of the particles along with tight packing specifications reduces losses in efficiency and performance due to band broadening.

Aeris Core-Shell Particles



Fully Porous Particles




The result is


- **3.6 µm core-shell particles** that can perform like sub-2 µm columns on both HPLC and UHPLC systems at a fraction of the pressure
- **1.7 µm core-shell particles** that can provide higher peak capacities compared to fully porous sub-2 µm columns on UHPLC systems

Selecting the Optimal Aeris Column for Your Applications


Aeris™ core-shell columns are designed for the separation of complex protein and peptide mixtures. Chromatographers can easily narrow down the column(s) that has a high probability of success for their separation by selecting from a variety of phase, pore size, and particle size options.



3.6 μm Core-Shell Particle



1.7 μm Core-Shell Particle




Aeris PEPTIDE

Recommended for the separation of low molecular weight peptides and peptide mapping.

- XB-C18 chemistry best suited for resolving peptides
- 1.7 μm and 3.6 μm particles for method development flexibility
- Small pore optimized for peptide diffusion


3.6 μm Core-Shell Particle



Aeris WIDEPORE

Recommended for the separation of intact proteins and large oligonucleotides.

- XB-C18, XB-C8, and C4 phases for alternate selectivities
- 3.6 μm particle for system flexibility
- Thin shell optimized for fast protein adsorption/desorption
- High pore permeability for improved separation of very large proteins (up to 400 kDa)




**Aeris WIDEPORE
XB-C18 and
Aeris PEPTIDE XB-C18
make a perfect pair
for peptide mapping.
See p. 32 for more details.**

Improve Resolution on ANY System by Leveraging Low Backpressure


Aeris™ PEPTIDE and Aeris WIDEPORE 3.6µm columns can **perform like sub-2µm columns at a fraction of the backpressure**. This allows chromatographers to utilize the resolving power of longer length (or coupled) columns without exceeding the pressure limits of their HPLC system. Scientists analyzing proteins and peptides can now have ultra-high resolution on HPLC or UHPLC systems.

Sub-2µm Performance at a Fraction of the Backpressure

Aeris WIDEPORE 3.6µm XB-C18



*Waters® ACQUITY® BEH300 1.7µm C18



Conditions for both columns:

Column: Aeris WIDEPORE 3.6µm XB-C18
ACQUITY® BEH300 1.7µm C18
Dimensions: 150 x 2.1 mm
Mobile Phase: A: Water with 0.1% TFA
B: Acetonitrile with 0.1% TFA
Gradient: A/B (65:35) for 3 min to A/B (35:65) over 30 min

Flow Rate: 0.3 mL/min
Temperature: 40 °C
Injection Volume: 20 µL
Instrument: Agilent® 1200SL
Detection: UV @ 214 nm (ambient)
Sample: BSA (Bovine Serum Albumin) Tryptic Digest

* ACQUITY and Waters are registered trademarks, and BEH Technology is a trademark of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.


Using a UHPLC system?

Try Aeris PEPTIDE 1.7µm columns for ultra-high efficiency peptide maps and stability up to 1,000 bar.


See page 11!



Utilize Long Columns to Maximize Resolution on UHPLC Systems




150 x 2.1 mm



150 x 2.1 mm

Moving to a
Longer Column

250 x 2.1 mm





Conditions for all columns:
Column: Aeris WIDEPORE 3.6 μ m XB-C18
 ACQUITY® BEH300 1.7 μ m C18
Dimensions: as noted in chromatogram
Mobile Phase: A: Water with 0.1 % TFA
 B: Acetonitrile
Gradient: A/B (90:10) for 5 min to A/B (50:50) over 45 min
Flow Rate: 0.2 mL/min
Temperature: 22°C
Injection Volume: 20 μ L
Instrument: Agilent® 1200SL
Detection: UV @ 210 nm
Sample: Degraded Myoglobin

* Waters and ACQUITY are registered trademarks, and BEH Technology is a trademark of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.

Achieve UHPLC Performance on HPLC Systems by Replacing 3 μm and 5 μm Columns

The innovative structure of 3.6 μm Aeris™ core-shell particles was specially designed to provide sub-2 μm performance at backpressures similar to fully porous 3 μm and 5 μm particles. Aeris columns can deliver increased resolution for existing protein and peptide separations performed on fully porous 3 μm and 5 μm columns, using the same HPLC system!


Upgrade Existing Methods on 3 μm and 5 μm Fully Porous Columns to Aeris Core-Shell Technology



Conditions for both columns:
Column: Aeris WIDEPORE 3.6 μm XB-C18
ZORBAX® 300SB 3.5 μm C18
Dimensions: 150 x 4.6 mm
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 3 min to A/B (35:65) over 30 min
Flow Rate: 1.5 mL/min
Temperature: 40 °C
Injection Volume: 20 μL
Instrument: Agilent® 1200SL
Detection: UV @ 210 nm (ambient)
Sample: Degraded Myoglobin

* Agilent and ZORBAX are registered trademarks of Agilent Technologies, Inc. Phenomenex is not affiliated with Agilent Technologies, Inc. Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.

**Improving your
current method is fast and
easy with an
Aeris core-shell column.**




Increase Efficiency on UHPLC Systems with Sub-2 μm Core-Shell Particles

For labs that have adopted higher pressure capable UHPLC instruments, AERIS PEPTIDE 1.7 μm core-shell columns are an excellent solution for ultra-high resolution peptide and peptide mapping separations. Core-shell particle technology combined with a sub-2 μm particle size results in extremely high efficiencies that scientists can use to pull apart critical peaks.

Ultra-High Resolution Achieved with 1.7 μm Core-Shell Technology

Conditions for both columns:
Column: AERIS PEPTIDE 3.6 μm XB-C18
 AERIS PEPTIDE 1.7 μm XB-C18
Dimensions: 150 x 2.1 mm
Part Nos.: 00F-4057-AN
 00F-4056-AN
Mobile Phase: A: Water with 0.1 % TFA
 B: Acetonitrile with 0.08 % TFA
Gradient: A/B (97:3) for 1.5 min to A/B (60:40) over 11 min to
 A/B (5/95) for 1 min
Flow Rate: 0.5 mL/min
Temperature: 40 °C
Injection Volume: 5 μL
Instrument: Agilent® 1200SL
Detection: UV @ 214 nm (ambient)
Sample: Alpha-Casein Tryptic Digest




Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.

Minimize Adsorption and Maximize Recoveries for Accurate Results

Aeris™ phase chemistries and bonding technology create a highly inert surface, leading to greatly reduced irreversible adsorption, higher recoveries, and sharper, narrower peaks, providing high quality and accurate results for each consecutive analysis.

Maximize Recoveries of Hydrophobic Proteins



Conditions for both columns:

Column: Aeris WIDEPORE 3.6 μm C4
ZORBAX® 300SB 3.5 μm C3
Dimensions: 150 x 2.1 mm
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) to A/B (35:65) over 45 min

Flow Rate: 0.3 mL/min
Temperature: 40 °C
Injection Volume: 20 μL
Instrument: Agilent® 1200
Detection: UV @ 214 nm (ambient)
Sample: Human Epidermal Growth Factor


* Agilent and ZORBAX are registered trademarks of Agilent Technologies, Inc. Phenomenex is not affiliated with Agilent Technologies, Inc. Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.

Tightly Controlled Quality for Reproducible Data



Every Aeris column and batch of media undergoes quality assurance tests for particle size distribution (both solid core and final particle), surface coverage, carbon load, pore diameter, pore size distribution, and other parameters to ensure **exceptional reproducibility for worry-free methods and confident results.**


Batch-to-Batch Reproducibility



App ID 19910

Column: Aeris PEPTIDE 1.7 μ m XB-C18
Dimensions: 50 x 4.6 mm
Part No.: 00B-4506-E0
Mobile Phase: A: Water with 0.1% Formic Acid
 B: Acetonitrile with 0.1% Formic Acid
Gradient: A/B (95:5) to A/B (5:95) over 4 min
Flow Rate: 1.85 mL/min
Temperature: 30 °C
Injection Volume: 0.4 μ L
Detection: UV @ 254 nm (ambient)
Sample: Selectivity Test Mixture

Batch-to-Batch Reproducibility



App ID 19908

Column: Aeris WIDEPOR 3.6 μ m XB-C18
Dimensions: 150 x 4.6 mm
Part No.: 00F-4482-E0
Mobile Phase: A: Water with 0.1% Formic Acid
 B: Acetonitrile with 0.085% Formic Acid
Gradient: A/B (95:5) to A/B (5:95) over 20 min
Flow Rate: 1.0 mL/min
Temperature: 40 °C
Injection Volume: 0.2 μ L
Detection: UV @ 210 nm (ambient)
Sample: Mouse IgG



Greater Method Flexibility with Specialty Surface Chemistries

Aeris™ WIDEPORE columns are available in three surface chemistries (XB-C18, XB-C8, C4) to satisfy applications of all types, ranging from sticky, intact proteins to complex protein digests.

Aeris PEPTIDE columns utilize the XB-C18 chemistry, as it is optimal for peptides and peptide mapping applications.

The unique, sterically protected XB surface ligands are designed by bonding bulky isobutyl chains beside the alkyl chains, and then fully end-capping the surface to cover any remaining exposed silanols.

An added benefit of XB chemistry is its high temperature stability, which allows one to use elevated column temperatures up to 90 °C for improved peak shape and recovery.




The Aeris WIDEPORE C4 phase does not use the XB chemistry, as shorter chain alkyl phases have higher bonding densities, thus providing steric hindrance. This means that chemical stability, inertness, and low bleed are maintained. The Aeris WIDEPORE C4 phase is an excellent complement to the other phases, and is also temperature stable to 90 °C



Long Column Lifetimes Under Extreme Method Conditions

Aeris columns provide temperature stability up to 90 °C, and pH stability from 1.5 - 9, giving ample flexibility for method development and excellent column lifetime.

Over 1,000 Injections at 90 °C




Column: Aeris WIDEPORE 3.6 µm XB-C18
Dimensions: 50 x 4.6 mm
Part No.: 00B-4282-E0
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 3 min, then to A/B (35:65) over 20 min

Flow Rate: 1.5 mL/min
Temperature: 90 °C
Injection Volume: 10 µL
Detection: UV @ 214 nm (ambient)
Sample: Apomyoglobin Digest

Low Column Bleed for Amplified Mass Spec (MS) Sensitivity

Aeris columns show no significant phase bleed under LC/MS conditions, making them very suitable for protein and peptide analysis. Chemists can be assured accurate, dependable, and consistent results, time and time again.

Virtually No LC/MS Bleed



Column: Aeris WIDEPORE 3.6 µm XB-C18
Dimensions: 50 x 2.1 mm
Part No.: 00B-4282-AN
Mobile Phase: A: Water with 0.1 % Formic Acid
B: Acetonitrile with 0.1 % Formic Acid
Gradient: A/B (95:5) for 2.5 min, to A/B (5:95) hold for 0.5 min, then re-equilibrate

Flow Rate: 0.5 mL/min
Temperature: 25 °C
Detection: MS (API 4000™)
Positive Ion Mode
Q1 scan from 75 to 800 amu
Sample: Blank

Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.

Aeris™ WIDEPORE Columns

for Intact Protein and Large Polypeptide Separations

Aeris WIDEPORE columns are packed with 3.6µm core-shell particles that are specially engineered with a thin porous shell, large pores, and sterically protected XB surface chemistry to address the inherent separation challenges of proteins and large peptides. This unique mix of features results in low backpressures, fast rates of diffusion, and excellent selectivity, generating exceptional chromatographic resolution on both HPLC and UHPLC systems.

Recommended for...

- Protein structural characterization
- Stability indicating assays
- Post-translational modification identification
- PEGylated proteins, antibodies, biosimilars, etc.
- Impurity profiling
- Alternate peptide map selectivity
- Large oligonucleotides






Table of Contents

Aeris WIDEPORE

p. 18 Easy Method Development with Three Selectivities

p. 20 Maximize HPLC and UHPLC Resolving Power


p. 22 Applications

“ The Aeris WIDEPORE column has given our company the opportunity to separate 2 forms of a protein (PEGylated & non-PEGylated). Prior to using Aeris the 2 peaks demonstrated little or no resolution. However by using the Aeris column the 2 peaks are separated by 5 minutes which is excellent. ”

-LARGE PHARMACEUTICAL COMPANY


“ Started using the Aeris WIDEPORE XB-C18 and XB-C8 for oligonucleotides and aptamers with excellent results! Very good peak shapes and excellent plate counts on these columns. Really nice to see all of the peaks present in the samples w/o a very long run time. Columns seem to be very stable and have very reasonable backpressures! ”

-HEALTHCARE PRODUCTS COMPANY




Easy Method Development with Three Selectivities

Aeris™ WIDEPORÉ 3.6µm Core-Shell Stationary Phases:



XB-C18
**Maximum hydrophobicity
recommended for:**

- Proteins
- Hydrophilic proteins
- PEGylated proteins
- High temperature separations
- Alternative selectivity for peptide mapping



XB-C8
**Moderate hydrophobicity
recommended for:**

- Proteins
- Moderately hydrophobic proteins
- Monoclonal antibodies
- Glycosylated proteins
- High temperature separations




C4
**Low hydrophobicity
recommended for:**


- Very large proteins
- Very hydrophobic proteins
- Membrane proteins
- Least retentive

Want more information on
the novel XB chemistry?


See page 14!



Easy Method Development with Three Selectivities



Because optimal separation conditions are different for each protein, we suggest evaluating all three surface chemistries to uncover the most suitable one for your separation. Once a phase is selected, the method can be further optimized with tweaks to the mobile phase, flow rate, gradient, or column dimension (length, internal diameter).









Conditions for all columns:

Column: Aeris WIDEPORE 3.6 μm (as noted)
Dimensions: 100 x 4.6 mm
Mobile Phase: A: Water with 0.1% TFA
 B: Acetonitrile with 0.1% TFA
Gradient: A/B (97:3) for 3 min to A/B (35:65) over 20 min
Flow Rate: 1.5 mL/min

Temperature: 40°C
Injection Volume: 10 μL
Instrument: Agilent® 1200
Detection: UV @ 214 nm (ambient)
Sample: Lysozyme (1 mg/mL)

Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.

Aeris Phase Selectivity Differences



Maximize HPLC and UHPLC Resolving Power with Unique 3.6 μm Core-Shell Particle

3.6 μm core-shell technology combined with inert surface chemistries and tight packing specifications results in Aeris™ WIDEPORE columns **delivering exceptional resolving power at significantly lower backpressures**. Chromatographers now have the ability to generate higher quality data than typically produced by columns packed with fully porous particles for every protein analysis – on HPLC or UHPLC systems.

Conditions for both columns:

Column: ACQUITY® BEH300 1.7 μm C4
Aeris WIDEPORE 3.6 μm C4

Dimensions: 150 x 2.1 mm

Mobile Phase: A: Water with 0.1% TFA
B: Acetonitrile with 0.1% TFA

Gradient: A/B (97:3) to A/B (35:65) over 45 min

Flow Rate: 0.3 mL/min

Temperature: 40 °C



Injection Volume: 10 μL

Instrument: Agilent® 1200

Detection: UV @ 214 nm (ambient)

Sample: Human Epidermal Growth Factor (EGF)


Performance Equivalent to sub-2 μm Particle at Low Backpressure



* Waters and ACQUITY are registered trademarks, and BEH Technology is a trademark of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.


Increase Column Length to Improve Resolving Power

150 x 2.1 mm




Aeris WIDEPORE

3.6 μ m XB-C18




150 mm Length
Zoom-In

250 x 2.1 mm



Aeris WIDEPORE

3.6 μ m XB-C18



250 mm Length
Zoom-In

Conditions for both columns:

Column: Aeris WIDEPORE 3.6 μ m XB-C18

Dimensions: as noted

Mobile Phase: A: Water with 0.1 % TFA

B: Acetonitrile with 0.1 % TFA

Gradient: A/B (97:3) for 3 min to A/B (35:65) over 30 min

Flow Rate: 0.3 mL/min



Temperature: 40 °C

Injection Volume: 25 μ L

Instrument: Agilent® 1200

Detection: UV @ 214 nm (ambient)

Sample: BSA (Bovine Serum Albumin) Digest




**Longer Column =
Increased Resolution**

Applications

Intact Protein Characterization


Biosimilar Impurity Quantitation



Column: Aeris™ WIDEPORE 3.6µm XB-C8
Dimensions: 150 x 4.6 mm
Part No.: 00F-4481-E0
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (70:30) to A/B (35:65) over 30 min

Flow Rate: 1.0 mL/min
Temperature: 22°C
Injection Volume: 5 µL
Instrument: Agilent® 1200
Detection: UV @ 214 nm (ambient)
Sample: Interferon alpha-2a

Protein Reduction



Column: Aeris WIDEPORE 3.6µm C4
Dimensions: 150 x 4.6 mm
Part No.: 00F-4486-E0
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 3 min to A/B (35:65) over 30 min
Flow Rate: 1.2 mL/min
Temperature: 22°C
Injection Volume: 20 µL
Instrument: Agilent 1200 SL
Detection: UV @ 214 nm (ambient)
Sample: RNase subject to reduction
100 % intact
20 % reduced
40 % reduced
60 % reduced
100 % reduced

Aeris WIDEPORE 3.6µm C4
successfully monitors peak
shifts due to differences in
protein shape




Applications


Intact Monoclonal Antibody (mAb) Separation



Human mAb



App ID 19846



App ID 19846


Column: Aeris WIDEPORE 3.6 μm XB-C18
Dimensions: 50 x 2.1 mm
Part No.: 00B-4482-AN

Mobile Phase: A: Water with 0.1 % Formic Acid
 B: Acetonitrile with 0.1 % Formic Acid
Gradient: A/B (90:10) to A/B (10:90) over 6 min

Step No.	Time(min)	% A	% B
1	0	90	10
2	0.7	66	34
3	5	55	45
4	6	10	90

Flow Rate: 0.5 mL/min
Temperature: 22 °C
Detection: UV @ 214 (ambient)
Sample: Monoclonal antibody

Clipped Variants




App ID 19845

Column: Aeris WIDEPORE 3.6 μm XB-C18
Dimensions: 250 x 4.6 mm
Part No.: 00G-4482-E0
Mobile Phase: A: Water with 0.1 % TFA
 B: Acetonitrile/IPA (50:50) with 0.1 % TFA
Gradient: A/B (90:10) to A/B (35:65) over 60 min
Flow Rate: 1.0 mL/min
Temperature: 22 °C
Injection Volume: 25 μL
Instrument: Agilent® 1200
Detection: UV @ 214 nm (ambient)
Sample: Proprietary customer monoclonal antibody
 with clipped variants




Applications

Post-Translational Modification Analysis

Oxidation



Deamidation



Column: Aeris™ WIDEPORE 3.6µm XB-C18
Dimensions: 100 x 4.6 mm
Part No.: 00D-4482-E0
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 3 min to A/B (15:85) over 45 min

Flow Rate: 1.2 mL/min
Temperature: 22°C
Injection Volume: 50 µL
Instrument: Agilent® 1100
Detection: UV @ 214 nm (ambient)
Sample: Insulin oxidized using 3% hydrogen peroxide

Column: Aeris WIDEPORE 3.6µm XB-C18
Dimensions: 150 x 4.6 mm
Part No.: 00F-4482-E0
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.085 % TFA
Gradient: A/B (90:10) to A/B (35:65) over 10 min


Flow Rate: 1.2 mL/min
Temperature: 40°C
Injection Volume: 1 µL
Instrument: Agilent® 1100
Detection: UV @ 214 nm (ambient)
Sample: Proprietary intact insulin 6 kDa deamidated

Applications

Peptide Mapping




Bovine Serum Albumin Tryptic Map



App ID 19835

Apomyoglobin Digest



App ID 19844

Column: Aeris WIDEPORE 3.6 μ m XB-C18
Dimensions: 250 x 2.1 mm
Part No.: 00G-4282-AN
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 3 min to A/B (35:65) over 47 min
Flow Rate: 0.3 mL/min
Temperature: 40 °C
Injection Volume: 10 μ L
Instrument: Agilent® 1200SL
Detection: UV @ 214 nm (ambient)
Sample: BSA Tryptic Digest

Column: Aeris WIDEPORE 3.6 μ m XB-C18
Dimensions: 150 x 4.6 mm
Part No.: 00F-4282-E0
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 3 min to A/B (35:65) over 30 min
Flow Rate: 1.5 mL/min
Temperature: 22 °C
Injection Volume: 20 μ L
Instrument: Agilent® 1200
Detection: UV @ 214 nm
Sample: Apomyoglobin Digest

Aeris™ PEPTIDE Columns

for Peptide and Peptide Mapping Separations

Based on core-shell particle technology, Aeris PEPTIDE particles are designed with small pores, inert XB-C18 surface chemistry, and two different particle sizes (3.6 μm and 1.7 μm) to meet the resolution demands of chromatographers performing complex peptide and peptide map separations on HPLC and/or UHPLC systems.

Aeris PEPTIDE columns are built for the following:

- Synthetic peptide impurity analysis
- Peptide mapping
- Identifying protein modifications
 - Glycosylation
 - Substitution
 - Truncation
- Analyzing post-translational modifications
 - Deamidation
 - Oxidation
 - Deletions






Table of Contents


Aeris PEPTIDE

- p. 28 Select the Most Suitable Aeris PEPTIDE Column**
- p. 29 Maximum Performance on UHPLC Systems**
- p. 30 Ultra-High Resolving Power on HPLC and UHPLC Systems**
- p. 32 Bundle Aeris PEPTIDE with Aeris WIDEPORE for Detailed Peptide Maps**
- p. 33 Applications**

“Using the Aeris PEPTIDE column gave us the same resolution separation for cyanobacterial peptides that we could achieve using a smaller particle size column, but with far lower back-pressures. This will allow us to transfer the methods to lower pressure HPLC systems whilst retaining our separation.”

-LARGE PHARMACEUTICAL COMPANY

Select the Most Suitable Aeris™ PEPTIDE Column to Achieve Your Separation Goals



The family of Aeris PEPTIDE XB-C18 columns is designed to provide versatility for the development of peptide separation methods. Depending on your resolution, throughput goals, and pressure capabilities of your system, you can choose from two particle sizes with unique performance attributes, as well as several column lengths to select the most suitable column for seamless method development and excellent results.

- UHPLC system required
- HPLC or UHPLC compatible

Maximize Performance on UHPLC Systems with Aeris PEPTIDE 1.7 μm Technology




With pressure stability up to 1,000 bar and the high efficiencies brought about by core-shell particle technology, the sub-2 μm Aeris PEPTIDE column produces breakthrough chromatographic performance on UHPLC systems. Use Aeris PEPTIDE 1.7 μm columns to boost the performance of sub-2 μm fully porous peptide mapping methods.


Increase Peak Count with 1.7 μm Aeris Core-Shell Technology

Conditions for both columns:
Column: Aeris PEPTIDE 1.7 μm XB-C18
ACQUITY® BEH300 1.7 μm C18
Dimensions: 150 x 2.1 mm
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.08 % TFA
Gradient: A/B (97:3) for 1.5 min to A/B (60:40) over 11 min to
A/B (5:95) over 1 min
Flow Rate: 0.5 mL/min
Temperature: 40 °C
Injection Volume: 5 μL
Instrument: Agilent® 1200SL
Detection: UV @ 214 nm (ambient)
Sample: Alpha-Casein Tryptic Digest

Aeris PEPTIDE 1.7 μm XB-C18



*Waters® ACQUITY® BEH300 1.7 μm C18





* Waters and ACQUITY are registered trademarks, and BEH Technology is a trademark of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.

Ultra-High Resolving Power on HPLC and UHPLC Systems with Aeris PEPTIDE 3.6 μm Columns

The Aeris™ PEPTIDE 3.6 μm core shell column was designed with one purpose in mind: to maximize the separation of large numbers of peptides on any HPLC or UHPLC system. Because core-shell particles remove the backpressure constraints of HPLC or UHPLC systems, chromatographers can **achieve the ultra-high performance of similar length sub-2 μm columns at a fraction of the backpressure.**

UHPLC Performance at HPLC Compatible Backpressures


Conditions for both columns:
Column: Aeris PEPTIDE 3.6 μm XB-C18
ACQUITY® BEH300 1.7 μm C18
Dimensions: 150 x 2.1 mm
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 5 min to A/B (55:45) over 55 min
Flow Rate: 0.3 mL/min
Temperature: 40 °C
Instrument: Agilent® 1200SL
Detection: UV @ 214 nm (ambient)
Sample: Human IgG Tryptic Digest



* Waters and ACQUITY are registered trademarks, and BEH Technology is a trademark of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Study was performed using new columns and, to the extent possible, identical experimental conditions were applied. Comparative separations may not be representative of all applications.



**Use longer (or coupled)
3.6 μm columns on
UHPLC and HPLC
systems to resolve
critical peaks**

For applications like peptide separations and peptide mapping where resolution is the primary goal, the lower backpressure of Aeris PEPTIDE 3.6 μm core-shell columns allow one to use longer columns for higher resolving power resulting in increased separation of closely eluting peptides.




Conditions for both columns:
Column: Aeris PEPTIDE 3.6 μm XB-C18
Dimensions: as noted
Mobile Phase: A: Water with 0.1 % Formic Acid
 B: Acetonitrile with 0.1 % Formic Acid
Gradient: A/B (97:3) for 3 min to A/B (35:65) over 30 min
Flow Rate: 1.2 mL/min
Temperature: 40 °C
Injection Volume: 25 μL
Instrument: Agilent® 1200
Detection: UV @ 214 nm (ambient)
Sample: Bovine Serum Albumin (BSA) Tryptic Digest

Utilize Long Columns to Maximize Separation Power



150 x 4.6 mm

140 bar



250 x 4.6 mm

200 bar


Minimal Increase in Backpressure

Bundle Aeris™ PEPTIDE with Aeris WIDEPORE for Detailed Peptide Maps


Aeris PEPTIDE 3.6µm XB-C18 and Aeris WIDEPORE 3.6µm XB-C18 are a “must-have” pair for chromatographers who analyze complex peptide mixtures. Because each has a unique pore size and surface area, they exhibit different selectivity. Protein chemists can take advantage of this diversity to achieve the critical resolution of target peptides in various regions of the map, thus simplifying their method development.

Utilize Differences in Small and Large Pore Size Selectivity for Optimal Resolution

Aeris PEPTIDE 3.6µm XB-C18



Aeris WIDEPORE 3.6µm XB-C18





Conditions for both columns:
Column: Aeris PEPTIDE 3.6µm XB-C18
Aeris WIDEPORE 3.6µm XB-C18
Dimensions: 150 x 4.6 mm
Part Nos.: 00F-4507-E0
00F-4482-E0
Mobile Phase: A: Water with 0.1 % TFA
B: Acetonitrile with 0.1 % TFA
Gradient: A/B (97:3) for 3 min to A/B (35:65) over 30 min
Flow Rate: 1.5 mL/min
Temperature: 40°C
Instrument: Agilent® 1200SL
Detection: UV @ 214 nm (ambient)
Sample: Alpha-1 Acid Glycoprotein Tryptic Digest

Applications



Peptide Mapping on Core-Shell Technologies

Aeris PEPTIDE vs. Other Core-Shell Columns



Conditions same for all columns:

Columns: Aeris PEPTIDE 1.7 μm XB-C18
Poroshell® 2.7 μm SB-C18 120 Å
Ascentis® Express Peptide 2.7 μm C18

Dimensions: 150 x 2.1 mm

Mobile Phase: A: Water with 0.1 % Formic Acid
B: Acetonitrile with 0.08 % Formic Acid

Gradient: A/B (97:3) for 1.5 min to A/B (60:40)
over 11 min to A/B (5:95) over 1 min


Flow Rate: 0.5 mL/min

Temperature: 40 °C

Instrument: Agilent® 1200SL

Detection: UV @ 214 nm (ambient)

Sample: Alpha-Casein Tryptic Digest




* Agilent and Poroshell are registered trademarks of Agilent Technologies, Inc. Ascentis Express Peptide is a registered trademark of Sigma-Aldrich Biotechnology. Phenomenex is not affiliated with Agilent Technologies, Inc or Sigma-Aldrich Biotechnology. Comparative separations may not be representative of all applications.

Extend the Lifetime of your Aeris Core-Shell Columns with SecurityGuard ULTRA

The SecurityGuard ULTRA guard cartridge system protects Aeris core-shell columns from damaging chemical contaminants, protein adsorption, and microparticulates. This innovative and easy-to-use column protection system will not alter chromatography or contribute to extra dead volume and is pressure rated up to 20,000 psi for UHPLC systems.




With SecurityGuard ULTRA



Virtually no change in chromatography when using SecurityGuard ULTRA!

Without SecurityGuard ULTRA



Conditions same for both separations:

Columns: Aeris WIDEPORE 3.6 μm XB-C18

Dimensions: 150 x 4.6 mm

Mobile Phase: A: Water with 0.1% TFA

B: Acetonitrile with 0.085% TFA

Gradient: A/B (97:3) for 3 min to A/B (35:65) over 30 min

Flow Rate: 1.2 mL/min

Temperature: 40 °C


Instrument: Agilent® 1200

Detection: UV @ 214 nm (ambient)


Sample: 1. Intact RNase A

2. Reduced RNase A

With SecurityGuard ULTRA



Without SecurityGuard ULTRA



Contaminant and particle buildup

guarantee

If SecurityGuard ULTRA cartridge protection system does not perform as well or better than your current guard cartridge system of similar phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

Ordering information on page 35

Ordering Information



Aeris WIDEPORE 3.6 µm Minibore Columns (mm)

	50 x 2.1	100 x 2.1	150 x 2.1	250 x 2.1	SecurityGuard™ ULTRA Cartridges*
					3/pk
XB-C18	00B-4482-AN	00D-4482-AN	00F-4482-AN	00G-4482-AN	AJO-8783
XB-C8	00B-4481-AN	00D-4481-AN	00F-4481-AN	00G-4481-AN	AJO-8785
C4	00B-4486-AN	00D-4486-AN	00F-4486-AN	00G-4486-AN	AJO-8899

Aeris WIDEPORE 3.6 µm Analytical Columns (mm)

	100 x 4.6	150 x 4.6	250 x 4.6	SecurityGuard™ ULTRA Cartridges*
				3/pk
XB-C18	00D-4482-E0	00F-4482-E0	00G-4482-E0	AJO-8769
XB-C8	00D-4481-E0	00F-4481-E0	00G-4481-E0	AJO-8771
C4	00D-4486-E0	00F-4486-E0	00G-4486-E0	AJO-8901

Aeris PEPTIDE 1.7 µm Minibore Columns (mm)

	50 x 2.1	100 x 2.1	150 x 2.1	SecurityGuard™ ULTRA Cartridges*
				3/pk
XB-C18	00B-4506-AN	00D-4506-AN	00F-4506-AN	AJO-8948

Aeris PEPTIDE 3.6 µm Minibore Columns (mm)

	50 x 2.1	100 x 2.1	150 x 2.1	250 x 2.1	SecurityGuard™ ULTRA Cartridges*
					3/pk
XB-C18	00B-4507-AN	00D-4507-AN	00F-4507-AN	00G-4507-AN	AJO-8948

Aeris PEPTIDE 3.6 µm Analytical Columns (mm)

	100 x 4.6	150 x 4.6	250 x 4.6	SecurityGuard™ ULTRA Cartridges*
				3/pk
XB-C18	00D-4507-E0	00F-4507-E0	00G-4507-E0	AJO-8946

* SecurityGuard ULTRA cartridges require holder part number, AJO-9000

SecurityGuard™ ULTRA Cartridge Holder* (for 2.1 to 4.6 mm ID columns)

SecurityGuard ULTRA Guard Cartridge Holder	ea
	AJO-9000

Material Characteristics

Packing Material	Total Particle Size (µm)	Porous Shell (µm)	Core Size (µm)	Pore Size (Å)	pH Stability	Temp Stability	Pressure Stability
Aeris WIDEPORE	3.6	0.2	3.2	200	1.5 - 9	90 °C	600 bar
Aeris PEPTIDE	1.7	0.22	1.25	100	1.5 - 9	90 °C	1000 bar
Aeris PEPTIDE	3.6	0.5	2.6	100	1.5 - 9	90 °C	600 bar



guarantee

If Aeris core-shell columns do not provide at least an equivalent separation as compared to a competing column of the same phase, return the column with comparative data within 45 days for a FULL REFUND.

Australia

t: 02-9428-6444
 f: 02-9428-6445
 auinfo@phenomenex.com

Austria

t: 01-319-1301
 f: 01-319-1300
 anfrage@phenomenex.com

Belgium

t: 02 503 4015 (French)
 t: 02 511 8666 (Dutch)
 f: +31 (0)30-2383749
 beinfo@phenomenex.com

Canada

t: (800) 543-3681
 f: (310) 328-7768
 info@phenomenex.com

Denmark

t: 4824 8048
 f: +45 4810 6265
 nordicinfo@phenomenex.com

Finland

t: 09 4789 0063
 f: +45 4810 6265
 nordicinfo@phenomenex.com

France

t: 01 30 09 21 10
 f: 01 30 09 21 11
 franceinfo@phenomenex.com

Germany

t: 06021-58830-0
 f: 06021-58830-11
 anfrage@phenomenex.com

India

t: 040-3012 2400
 f: 040-3012 2411
 indiainfo@phenomenex.com

Ireland

t: 01 247 5405
 f: +44 1625-501796
 eireinfo@phenomenex.com

Italy

t: 051 6327511
 f: 051 6327555
 italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
 f: +31 (0)30-2383749
 nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
 f: 001-310-328-7768
 tecnicomx@phenomenex.com

The Netherlands

t: 030-2418700
 f: 030-2383749
 nlinfo@phenomenex.com

New Zealand

t: 09-4780951
 f: 09-4780952
 nzinfo@phenomenex.com

Norway

t: 810 02 005
 f: +45 4810 6265
 nordicinfo@phenomenex.com

Puerto Rico

t: (800) 541-HPLC
 f: (310) 328-7768
 info@phenomenex.com

Sweden

t: 08 611 6950
 f: +45 4810 6265
 nordicinfo@phenomenex.com

United Kingdom

t: 01625-501367
 f: 01625-501796
 ukinfo@phenomenex.com

United States

t: (310) 212-0555
 f: (310) 328-7768
 info@phenomenex.com

**All other countries:
Corporate Office USA** 

t: (310) 212-0555
 f: (310) 328-7768
 info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Aeris and SecurityGuard are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks, and BEH Technology is a trademark of Waters Corporation. Agilent, ZORBAX, and Poroshell are registered trademarks of Agilent Technologies, Inc. API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license. Ascentis is a registered trademark of Sigma-Aldrich Biotechnology LP and Sigma-Aldrich Co.

Disclaimer

Phenomenex is not affiliated with Agilent Technologies, Inc, Sigma-Aldrich Biotechnology LP and Sigma-Aldrich Co., or Waters Corporation. Comparative separations may not be representative of all applications.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

© 2012 Phenomenex, Inc. All rights reserved.

Core-Shell Technology
for Proteins and Peptides

