

Organic Size Exclusion/Gel Permeation for Polymer Analysis

- 5 and 10 μm particle sizes
- Narrow bore (4.6 mm ID) solvent-saver to preparative columns available
- Alternative to Agilent® (Polymer Labs) PLgel™, Waters® Styragel® and Ultrastyrigel™, and other columns (see p. 239)
- Highly cross-linked for mechanical and chemical stability
- Temperature stable to 140 °C

Phenogel is available in seven different pore sizes, ranging from 50 Å to 10⁶ Å†, and a linear bed configuration. Pore size distribution and pore volume are closely controlled parameters in the manufacturing process accounting for the high resolution, tight linear calibration curves, and excellent column-to-column reproducibility.


Sample Elution

Each standard dimension Phenogel column (300 x 7.8 mm) has an internal volume of 15 mL that is distributed as follows:

- 3 mL is occupied by the solid portions of the gel particles (20% of total column volume)
- 6 mL is the pore volume of the packing material (40% of total column volume)
- 6 mL is the interstitial volume or volume between the gel particles (40% of total column volume)

Thus, about 6 mL of solvent must elute through each column before even the largest molecules can emerge, while the smallest molecules emerge with the total column volume of 12 mL. This constant distribution of volume makes it possible to predict the amount of solvent and time necessary to complete any analysis.

SEM Photos of Phenogel Polymer Beads



If Phenogel analytical columns do not provide at least equivalent separation as compared to a competing column of the similar particle size, phase, and dimensions, send in your comparative data within 45 days and keep the column for FREE.

Technical Specifications

Material:	SDVB
Particle Size:	5, 10 μm
Porosities:	50 Å to 10 ⁶ Å†, and mixed beds
Maximum Pressure:	1500 psi
Maximum Temperature:	140 °C
Minimum Efficiency*:	5 μm : 45,000 p/m** 10 μm : 35,000 p/m**
Typical Flow Rates:	4.6 mm ID: 0.35 mL/min 7.8 mm ID: 1.0 mL/min 21.2 mm ID: 7.0 mL/min
End Fittings:	Valco® Compatible


* Tested in THF ** For 300 x 7.8 mm ID columns

† See note on p. 391 regarding pore sizes and exclusion limits

Column Selection by Molecular Weight

Sample Type	Molecular Weight	Phenogel Column
Small Organics	100 - 3 K	50 Å
	500 - 6 K	100 Å
	1 K - 15 K	500 Å
Resins	1 K - 75 K	10 ³ Å
	5 K - 500 K	10 ⁴ Å
	10 K - 1,000 K	10 ⁵ Å
High MW Polymers	60 K - 10,000 K	10 ⁶ Å
	100 - 10,000 K	Linear(2)

Column Molecular Weight Calibration Curves



Solvent and Temperature Compatibility

- Phenogel columns are packed in tetrahydrofuran (THF)
- Columns can also be shipped in solvents such as DMF, Methylene Chloride, NMP, and *o*-CP to help minimize equilibration time

Solvent Compatibility Table

Mobile Phase Solvent	Phenogel Pore Size:								Linear & Mixed	Suggested Operating Temp.
	50 (Å)	100	500	10 ³	10 ⁴	10 ⁵	10 ⁶			
Acetone	Y	Y	Y	Y	Y	Y	Y	Y		
Benzene	Y	Y	Y	Y	Y	Y	Y	Y		
Carbon Tetrachloride	Y	Y	Y	Y	Y	Y	Y	Y		
Chloroform	Y	Y	Y	Y	Y	Y	Y	Y		
30 % HFIP/Chloroform	Y	Y	Y	Y	Y	Y	Y	Y		
Diethyl Ether	Y	Y	Y	Y	Y	Y	Y	Y		
Dimethylacetamide (DMAC)	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C	
Dimethylformamide (DMF)	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C	
Dioxane	Y	Y	Y	Y	Y	Y	Y	Y		
DMSO	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C	
Ethyl Acetate	Y	Y	Y	Y	Y	Y	Y	Y		
Hexafluoroisopropanol (HFIP)	Y	Y	Y	Y	Y	Y	Y	Y		
Hexane	Y	Y	Y	Y	Y	Y	Y	Y		
M-Cresol	Y*	Y	Y	Y	Y	Y	Y	Y	100 °C	
Methyl Ethyl Ketone	Y	Y	Y	Y	Y	Y	Y	Y		
Methylene Chloride	Y	Y	Y	Y	Y	Y	Y	Y		
<i>o</i> -Chlorophenol	Y*	Y	Y	Y	Y	Y	Y	Y	100 °C	
<i>o</i> -Dichlorobenzene	Y*	Y	Y	Y	Y	Y	Y	Y	135 °C	
Quinolin	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C	
Tetrahydrofuran	Y	Y	Y	Y	Y	Y	Y	Y		
Toluene	Y	Y	Y	Y	Y	Y	Y	Y		
Trichlorobenzene	Y*	Y	Y	Y	Y	Y	Y	Y	135 °C	
Water	N	N	N	N	N	N	N	N		
Xylene	Y	Y	Y	Y	Y	Y	Y	Y		

*Not recommended on 5 μm 50 Å columns.

N = Not Compatible
Y = Compatible




Solvent Switching Considerations

Although Phenogel columns are rugged and can withstand strong solvent changes, care should be exercised when switching from high-swell solvents (A) to low-swell solvents (B, C, and D). Improper solvent switches can result in a void. Best results are attained when an intermediate-swell solvent is used and column lifetime is improved. Contact Phenomenex regarding solvents not listed below.

Column life can be maximized by dedicating certain columns to certain solvents. This will also minimize solvent switches. If care is not taken, a void may occur.

- Reduce flow rate to 0.2 mL/min
- Backpressure must NEVER exceed 1500 psi
- Always check solvent miscibility in a beaker or follow the solvent miscibility table on page 389 before proceeding with ANY solvent switch.
- Compare the swell characteristics of solvent 1 (old solvent) to solvent 2 (new solvent) and use the following guidelines:
 - If solvent 1 and solvent 2 belong to the same swell category (see table below), check the solvent miscibility and proceed with the switch.
 - If solvent 1 and solvent 2 belong to successive swell categories as indicated by the arrows in the table below, check the miscibility and proceed with the switch.
 - If solvent 1 and solvent 2 DO NOT belong to the same OR successive swell categories, switch to an intermediate solvent FIRST, as indicated by the arrows in the table.




Phenogel™ SEC/GPC Columns

Phenogel and Excipients Analysis

Gel permeation chromatography using Phenogel columns is an excellent method for measuring the molecular weight distribution and lot-to-lot consistency of fillers and dispersants.


Polyethylene Glycol 330

Column: Phenogel 5 μm 50 Å, 100 Å, 500 Å
Dimensions: 300 x 7.8 mm
Solvent: THF
Flow Rate: 1.0 mL/min
Detection: Differential Refractometer
Injection Volume: 100 μL 0.25% w/v
Temperature: Ambient
Sample: 1. dp7 546 MW 5. dp3 194 MW
2. dp6 458 MW 6. dp2 106 MW
3. dp5 370 MW 7. dp1 62 MW
4. dp4 282 MW




Polyethylene Glycol 106

Column: Phenogel 5 μm 50 Å, 100 Å, 500 Å
Dimensions: 300 x 7.8 mm
Mobile Phase: THF
Flow Rate: 1 mL/min
Detection: Refraction Index (RI) (40 °C)
Temperature: 40 °C
Sample: 1. PEG 106
2. API peak A (unknown)
3. API peak B (unknown)



Polyvinylpyrrolidone

Column: Phenogel 5 μm Linear(2) x2
Dimensions: 300 x 7.8 mm
Part No: 00H-3259-KO
Mobile Phase: 10 mM Lithium bromide in DMF
Flow Rate: 2 mL/min
Detection: Refraction Index (RI) (40 °C)
Column Temp: 40 °C
Sample: 1. Polyvinylpyrrolidone (PVP)




Phenogel™ SEC/GPC Columns

50 Å - 10⁶ Å Columns

- High resolution at low cost
- Customize your analysis by coupling different pore-size columns
- Wide range of solvent compatibility


Polymethyl Methacrylates (Wide MW Range)

Column: Phenogel 5 μm 10⁶ Å, 10⁴ Å, 10³ Å, 500 Å
Dimensions: 300 x 7.8 mm
Solvent: THF
Flow Rate: 1.0 mL/min
Detection: Differential Refractometer
Injection Volume: 100 μL 0.25% w/v
Sample:




Closely Related Hydrocarbons

Column: Phenogel 5 μm 50 Å, 100 Å, 500 Å
Dimensions: 300 x 7.8 mm
Solvent: THF
Flow Rate: 1.0 mL/min
Detection: Differential Refractometer
Injection Volume: 100 μL 0.25% w/v
Temperature: Ambient
Sample:




Polyethylene Oxide (PEO)

Column: Phenogel 10 μm 10⁵, 10⁴, 10³ Å
Dimensions: 300 x 7.8 mm
Mobile Phase: DMF (0.1 M LiBr)
Flow Rate: 1.0 mL/min
Detection: Differential Refractometer
Injection Volume: 100 μL 0.125% w/v
Temperature: 50 °C
Sample: 400,000 MW




Polystyrenes (Wide MW Range)

Column: Phenogel 10 μm 10⁵, 10⁴, 10³ Å
Dimensions: 300 x 7.8 mm
Mobile Phase: THF
Flow Rate: 1.0 mL/min
Detection: Differential Refractometer
Injection Volume: 100 μL 0.125% w/v
Temperature: Ambient
Sample:




Poly-(α-Methyl Styrene) (Wide MW Range)

Column: Phenogel 5 μm 10⁵, 10⁴, 10³, 500 Å
Dimensions: 300 x 7.8 mm
Solvent: THF
Flow Rate: 1.0 mL/min
Detection: Differential Refractometer
Injection Volume: 100 μL 0.25% w/v
Sample:



Polybutadienes (Wide MW Range)

Column: Phenogel 5 μm 10⁵, 10⁴, 10³, 500 Å
Dimensions: 300 x 7.8 mm
Solvent: THF
Flow Rate: 1.0 mL/min
Detection: Differential Refractometer
Injection Volume: 100 μL 0.25% w/v
Sample:



Phenogel™ SEC/GPC Columns


Linear Columns

- Linear calibration to 10 million daltons
- Long column lifetime
- Excellent mechanical stability
- Excellent for analyzing a wide range of molecular weights




Mixed Polystyrene Standard

Column: Phenogel 5 µm Linear(2)
Dimensions: 300 x 7.8 mm
Part No.: 00H-3259-K0
Mobile Phase: THF
Flow Rate: 1.0 mL/min
Detection: RI
Injection Volume: 50 µL
Temperature: 35 °C
Sample: Polystyrene standards injected
 1. 2,860,000 MW
 2. 277,000 MW
 3. 9,350 MW



Calibration Curve: Linear (2) - Phenogel 5 µm 300 x 7.8 mm



Narrow Bore Columns


An Improved Dimension in GPC Analysis

- Decrease solvent consumption
- Retain same elution profile
- Reduce solvent disposal costs

Phenogel-NB (Narrow Bore) columns are optimized to reduce solvent consumption. The Phenogel-NB columns have a 4.6 mm column ID and run at 0.35 mL/min, reducing solvent consumption and disposal costs up to 65 %!

Loading

With narrow bore GPC/SEC columns, the volume in which the sample elutes is significantly decreased, thus increasing the effective concentration of the sample. This increase in sensitivity is exploited in HPLC, but in GPC it leads to overloading effects and proportionally lower sample loadings must be used.



Phenogel™ SEC/GPC Columns

guarantee

If Phenogel analytical columns do not provide at least equivalent separation as compared to a competing column of the similar particle size, phase, and dimensions, send in your comparative data within 45 days and keep the column for FREE.

Ordering Information

5 µm Columns (mm)		Guards	
		300 x 7.8	50 x 7.8
Pore Size	MW Range		
50 Å	100-3 K	00H-0441-KO	03B-2088-KO
100 Å	500-6 K	00H-0442-KO	03B-2088-KO
500 Å	1 K-15 K	00H-0443-KO	03B-2088-KO
10 ³ Å	1 K-75 K	00H-0444-KO	03B-2088-KO
10 ⁴ Å	5 K-500 K	00H-0445-KO	03B-2088-KO
10 ⁵ Å	10 K-1,000 K	00H-0446-KO	03B-2088-KO
10 ⁶ Å	60 K-10,000 K	00H-0447-KO	03B-2088-KO
		300 x 7.8	50 x 7.8
Mixed Beds			
Linear(2)	100-10,000 K	00H-3259-KO	03B-2088-KO

5 µm Narrow Bore (NB) Columns (mm)		Guards	
		300 x 4.6	30 x 4.6
Pore Size	MW Range		
50 Å	100-3 K	00H-0441-E0	03A-2088-E0
100 Å	500-6 K	00H-0442-E0	03A-2088-E0
500 Å	1 K-15 K	00H-0443-E0	03A-2088-E0
10 ³ Å	1 K-75 K	00H-0444-E0	03A-2088-E0
10 ⁴ Å	5 K-500 K	00H-0445-E0	03A-2088-E0
		300 x 4.6	30 x 4.6
Mixed Beds			
Linear(2)	100-10,000 K	00H-3259-E0	03A-2088-E0

10 µm Columns (mm)		Guards	
		300 x 7.8	50 x 7.8
Pore Size	MW Range		
50 Å	100-3 K	00H-0641-KO	03B-2090-KO
100 Å	500-6 K	00H-0642-KO	03B-2090-KO
500 Å	1 K-15 K	00H-0643-KO	03B-2090-KO
10 ³ Å	1 K-75 K	00H-0644-KO	03B-2090-KO
10 ⁴ Å	5 K-500 K	00H-0645-KO	03B-2090-KO
10 ⁵ Å	10 K-1,000 K	00H-0646-KO	03B-2090-KO
10 ⁶ Å	60 K-10,000 K	00H-0647-KO	03B-2090-KO
		300 x 7.8	50 x 7.8
Mixed Beds			
Linear(2)	100-10,000 K	00H-3260-KO	03B-2090-KO



All other column dimensions available. Phenogel columns are routinely shipped in THF. However, columns are also available in commonly used solvents such as Toluene and Chloroform as well as DMF, NMP, and other solvents. Refer to the chart below for the additional charge for these shipping solvents. Please specify shipping solvent when ordering.

Other Shipping Solvents:	
Methanol, Methylene Chloride, Cyclohexane, Ethyl Acetate, NMP, DMAC, DMF	
Size (mm)	Price
30 x 4.6	
50 x 7.8	
300 x 4.6	
300 x 7.8	

NOTE: Phenogel columns are routinely shipped in THF. Columns can be shipped in Toluene and Chloroform upon request at no additional charge.

Recommended Alternatives

Manufacturer	Columns
Agilent® (Polymer Labs)	PLgel™
Jordi Labs	Jordi Gel™ DVB Jordi Gel DVB Fluorinated Jordi Gel DVB Glucose
Polymer Standards Service (PSS)	SDV GRAM PolarSil PFG POLEFIN
Shodex®	GPC K-800 Series GPC KF-800 Series GPC KD-800 Series
Tosoh Bioscience®	TSKgel® SuperMultiporeHZ TSKgel SuperHZ TSKgel Hxl TSKgel SuperH TSKgel Hhr
Waters®	Styragel® UltraStyragel™ ACQUITY® APC™



For Column Heaters, see p. 366