

Clarity® BioSolutions for Synthetic DNA/RNA

Clarity Oligo-RP™ HPLC Columns

Reversed Phase LC for Purification and Characterization

- Easily separate N-1 failure sequences from target oligo with > 90 % purities
- Trityl-off purification of DNA, RNA, Thioates, and modified/labeled oligonucleotides
- Preparative dimensions and particle sizes for loads > 5 μ mole
- Purify oligos up to 60 mer in length
- Excellent column for reversed phase HPLC quality control (QC) testing


Clarity Oligo-RP has been specifically designed for the reversed phase purification of oligonucleotides with balanced hydrophobicity and polar selectivity. The media is based on composite particle TWIN™ technology. This technology gives improved selectivity and efficiency for oligonucleotides when compared to other hybrid, polymer, and silica particles found in the marketplace. It is available in 3, 5, and 10 μ m particle sized beads and in a variety of dimensions.

Preparative Purification on Oligo-RP


Reversed phase separation of oligonucleotides has advantages over other modes of separations such as ion-exchange. The Oligo-RP phase allows high loadability and delivers high recovery and purity, eliminating the need for extra purification steps. This is achieved through an ion-pair separation of the trityl-off oligonucleotide from failure products and other impurities.

DNA Purification (A) Preparative (B) Analytical QC

Column: Clarity 3 μ m Oligo-RP C18
Dimensions: (A) 50 x 10.0 mm
(B) 50 x 4.6 mm
Part No.: (A) 00B-4441-NO
(B) 00B-4441-EO
Mobile Phase: A: 50 mM TEAA pH 7.5 / 5 % Acetonitrile
B: Methanol
Gradient: 10 % to 60 % B in 20 minutes
Flow Rate: (A) 4.7 mL / min
(B) 1.0 mL / min
Detection: UV @ 260 nm
Sample: 20nt DNA



A 200 μ g (1 μ mole) 20mer DNA sample was loaded onto a 10 mm ID Clarity Oligo-RP column. Impurities were separated from the target sequence.




A Clarity Oligo-RP analytical column was used to verify the purity of the preparative purification. A purity of 92 % with a yield of 85 % was determined.

Separate N-1 Failure Sequences from Target N Sequence

The Oligo-RP sorbent is specifically designed to accommodate all possible interactive features of nucleosides with matching modes of reactivity to its own. The sorbent possesses hydrophobic, dipolar, π - π , and hydrogen bond donor/acceptor sites; this combination of interaction along with an ion-pairing reagent elicits a high degree of differential selectivity between nucleic acids. Thus it can recognize even the slightest changes in nucleotide sequence, such as a difference of one base (N and N-1) or substitution of one base for another.

DNA Purification of Failure N-1 from Target N Sequence


Column: Clarity 3 μ m Oligo-RP C18
Dimensions: 50 x 4.6 mm
Part No.: 00B-4441-EO
Mobile Phase: A: 50 mM TEAA pH 7.5
B: Methanol
Gradient: 10 % to 45 % B in 30 minutes
Flow Rate: 1 mL / min
Detection: UV @ 260 nm
Sample: 1. 40nt DNA with sequence
CTTCTGAACAGTTGATCTATGCACCTTCAGACTTATGATCA (2.5 μ g)
2. 39nt DNA with sequence
TTCTGAACAGTTGATCTATGCACCTTCAGACTTATGATCA (2.5 μ g)



Clarity Oligo-RP successfully separates a 40mer from a 39mer DNA oligonucleotide due to its excellent efficiency and resolving power.

Fingerprint of 40mer DNA

Column: Clarity 3 μ m Oligo-RP C18
Dimensions: 50 x 4.6 mm
Part No.: 00B-4441-EO
Mobile Phase: A: 50 mM TEAA pH 7.5 / 5 % Acetonitrile
B: Methanol
Gradient: 20 % to 25 % B in 20 minutes; hold at 5 minutes @ 25 % B
Flow Rate: 1 mL / min
Detection: UV @ 260 nm
Sample: 40nt DNA with sequence
5'-CTC CTG GGC AGT GGA TCT GCG CACTTC AGG CTC CTG GGC A-3'



Due to the high efficiency of the sorbent and ion-pairing interactions, a fingerprint of a crude 40mer DNA on Clarity Oligo-RP is produced illustrating baseline resolution of impurities from the final product.