


## Примеры хроматограмм

### Phthalate Esters: EPA Method 606



**Column:** Zebtron ZB-5  
**Dimensions:** 30 meter x 0.53 mm x 1.50  $\mu$ m  
**Part No.:** 7HK-G002-28  
**Injection:** Split 20:1 @ 300 °C, 1  $\mu$ L  
**Carrier Gas:** Helium @ 12.9 mL/ min (constant flow)  
**Oven Program:** 40 °C for 6 min to 300 °C @ 10 °C/min for 15 min  
**Detector:** FID @ 300 °C  
**Sample:**


1. Dimethyl Phthalate
2. Diethyl Phthalate
3. Di-n-butyl Phthalate
4. Butyl Benzyl Phthalate
5. Bis(2-ethylhexyl) Phthalate
6. Di-n-octyl Phthalate

## Исследования

### Nonhalogenated Organics Using GC/FID: EPA Method 8015B

**Column:** Zebtron ZB-5  
**Dimensions:** 30 meter x 0.53 mm x 1.50  $\mu$ m  
**Part No.:** 7HK-G002-28  
**Injection:** Split 20:1 @ 225 °C, 1  $\mu$ L  
**Carrier Gas:** Helium @ 4.8 mL/min (constant flow)  
**Oven Program:** 45 °C for 3 min to 60 °C at 5 °C/min to 275 °C at 20 °C/min  
**Detector:** FID @ 350 °C  
**Sample:** All analytes at 100  $\mu$ g/mL

1. Methanol (solvent)
2. Ethanol
3. Acetonitrile
4. Diethyl ether
5. Propanenitrile
6. 2-Butanone [MEK]
7. 2-Methyl-2-propanenitrile
8. Isobutanol
9. 1,4-Dioxane
10. Methylmethacrylate
11. 4-Methyl-2-Pentanone [MIBK]
12. Impurity
13. Ethylmethacrylate
14. Acrylamide




## Продукты питания, ароматические и вкусовые добавки

### Vitamin E and Sterols

**Column:** Zebron ZB-5  
**Dimensions:** 30 meter x 0.25 mm x 0.10  $\mu$ m  
**Part No.:** 7HG-G002-02  
**Injection:** Splitless @ 220 °C, 1  $\mu$ L  
**Carrier Gas:** Helium @ 1.8 mL/min (constant flow)  
**Oven Program:** 110 °C for 0.2 min to 140 °C @ 30 °C/min to 230 °C @ 10 °C/min for 6 min to 340 °C @ 10 °C/min for 15.8 min  
**Detector:** FID @ 340 °C  
**Sample:**

1. Squalene	8. $\gamma$ -Tocomoenoel
2. FFA C24:0	9. Stigmasta-3,5,-diene
3. $\delta$ -Tocopherol	10. Cholesterol
4. $\delta$ -Tocomoenoel	11. $\alpha$ -Tocopherol
5. Campesta-3,5,-diene	12. Campesterol
6. $\gamma$ -Tocopherol	13. Stigmasterol
7. Stigmasta-3,5,22,-triene	14. $\beta$ -Sitosterol




## Химическая промышленность

### Polynuclear Aromatic Hydrocarbons (PAHs)

**Column:** Zebron ZB-5 with 5 meter Guardian™  
**Dimensions:** 30 meter x 0.25 mm x 0.25  $\mu$ m  
**Part No.:** 7HG-G002-11-GGA  
**Injection:** Split 15:1 @ 260 °C, 1  $\mu$ L  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 170 °C to 240 °C @ 15 °C/min to 275 °C @ 4 °C/min to 320 °C @ 10 °C/min for 5 min  
**Detector:** MSD @ 280 °C; 50-450 amu  
**Sample:** All analytes are 100  $\mu$ g/mL in methylene chloride


1. Naphthalene	10. Chrysene
2. Acenaphthylene	11. Benzo[b]fluoranthene
3. Acenaphthene	12. Benzo[k]fluoranthene
4. Fluorene	13. Benzo[a]pyrene
5. Phenanthrene	14. Indeno[1,2,3-cd]pyrene
6. Anthracene	15. Dibenz[a,h]anthracene
7. Fluoranthene	16. Benzo[g,h,i]perylene
8. Pyrene	
9. Benz[a]anthracene	



### Aliphatic Amines

**Column:** Zebron ZB-5  
**Dimensions:** 30 meter x 0.25 mm x 0.25  $\mu$ m  
**Part No.:** 7HG-G002-11  
**Injection:** Split 15:1 @ 220 °C, 1  $\mu$ L  
**Carrier Gas:** Helium @ 1.0 mL/min (constant flow)  
**Oven Program:** 40 °C for 3 min to 130 °C @ 30 °C/min to 220 °C @ 15 °C/min  
**Detector:** MSD; 50-450 amu  
**Sample:** 1.58 mg/mL each

1. tert-Butylamine	12. Dibutylamine
2. sec-Butylamine	13. Octylamine
3. Isobutylamine	14. Nonylamine
4. n-Butylamine	15. Dipentylamine
5. Diisopropylamine	16. Tributylamine
6. Triethylamine	17. Decylamine
7. Isopentylamine	18. Contaminant
8. Pentylamine	19. Dihexylamine
9. Hexylamine	20. Dicyclohexylamine
10. Cyclohexylamine	21. Contaminant
11. Heptylamine	



Note: Peak 17 is reduced in size because the molecular ion peak is outside the selected mass scan range.