

Strata® Traditional SPE


guarantee

If Strata SPE products do not perform as well or better than your current SPE product of similar phase, mass and size, send in your comparative data within 45 days and keep the Strata SPE product for FREE!

Mixed-Mode Sorbents

Screen-C

Incorporates the hydrophobic selectivity of a C8 phase and strong cation-exchange for the extraction of basic drugs from biological matrices.



Ordering Information

Format	Sorbent Mass	Part Number	Unit	Price
Tube				
	100 mg	8B-S016-EAK**	1 mL (100/box)	
	100 mg	8B-S016-EBJ	3 mL (50/box)	
	150 mg	8B-S016-SBJ	3 mL (50/box)	
	150 mg	8B-S016-SCH	6 mL (30/box)	
	200 mg	8B-S016-FBJ	3 mL (50/box)	
	200 mg	8B-S016-FCH	6 mL (30/box)	
	300 mg	8B-S016-RBJ	3 mL (50/box)	
	300 mg	8B-S016-RCH	6 mL (30/box)	
	500 mg	8B-S016-HCH	6 mL (30/box)	

Giga™ Tube

	2 g	8B-S016-KDG	12 mL (20/box)	
--	-----	-------------	----------------	--

96-Well Plate

	50 mg	8E-S016-DGB	2 Plates/Box	
	100 mg	8E-S016-EGB	2 Plates/Box	

Screen-C GF

Offers the selectivity of Screen-C in a gravity flow particle size for viscous samples.

Ordering Information

Format	Sorbent Mass	Part Number	Unit	Price
Tube				
	200 mg	8B-S026-FBJ	3 mL (50/box)	
	200 mg	8B-S026-FCH	6 mL (30/box)	
	500 mg	8B-S026-HBJ	3 mL (50/box)	
	500 mg	8B-S026-HCH	6 mL (30/box)	
	1 g	8B-S026-JCH	6 mL (30/box)	


Giga Tube

	500 mg	8B-S026-HDG	12 mL (20/box)	
	2 g	8B-S026-KEG	20 mL (20/box)	
	5 g	8B-S026-LFF	60 mL (16/box)	

**Tab-less tubes available. Contact Phenomenex for details.

ABW

Offers a strong cation-exchange group and a weak anion-exchange group for the extraction or fractionation of complex mixtures.




Ordering Information

Format	Sorbent Mass	Part Number	Unit	Price
Tube				
	200 mg	8B-S030-FBJ	3 mL (50/box)	
	1 g	8B-S030-JCH	6 mL (30/box)	
Giga Tube				
	2 g	8B-S030-KDG	12 mL (20/box)	
	5 g	8B-S030-LEG	20 mL (20/box)	
	10 g	8B-S030-MFF	60 mL (16/box)	
	20 g	8B-S030-VFF	60 mL (16/box)	

Screen-A

Incorporates the hydrophobic selectivity of a C8 phase and strong anion-exchange for the extraction of acidic drugs from biological matrices.



Ordering Information

Format	Sorbent Mass	Part Number	Unit	Price
Tube				
	100 mg	8B-S019-EAK	1 mL (100/box)	
	200 mg	8B-S019-FBJ	3 mL (50/box)	
	200 mg	8B-S019-FCH	6 mL (30/box)	
	500 mg	8B-S019-HCH	6 mL (30/box)	
96-Well Plate				
	25 mg	8E-S019-CGB	2 Plates/Box	