

# Strata® Traditional SPE


guarantee

If Strata SPE products do not perform as well or better than your current SPE product of similar phase, mass and size, send in your comparative data within 45 days and keep the Strata SPE product for FREE!

## Normal Phase Sorbents


### Silica (Si-1)

Unbonded silica particle that offers strong polar selectivity.



### NH<sub>2</sub>/WAX

This amino phase offers strong polar selectivity and hydrogen bonding under normal phase conditions or can be used as a weak anion-exchange sorbent.



#### Ordering Information

Format	Sorbent Mass	Part Number	Unit	Price
<b>Tube</b>				
	100 mg	8B-S012-EAK	1 mL (100/box)	
	200 mg	8B-S012-FBJ	3 mL (50/box)	
	500 mg	8B-S012-HBJ**	3 mL (50/box)	
	500 mg	8B-S012-HCH**	6 mL (30/box)	
	1 g	8B-S012-JCH**	6 mL (30/box)	

#### Giga™ Tube

	500 mg	8B-S012-HDG	12 mL (20/box)	
	1 g	8B-S012-JDG	12 mL (20/box)	
	2 g	8B-S012-KDG	12 mL (20/box)	
	5 g	8B-S012-LEG	20 mL (20/box)	
	10 g	8B-S012-MFF	60 mL (16/box)	
	20 g	8B-S012-VFF	60 mL (16/box)	
	50 g	8B-S012-YSN	150 mL (8/box)	
	70 g	8B-S012-ZSN	150 mL (8/box)	

#### 96-Well Plate

	50 mg	8E-S012-DGB	2 Plates/Box	
	100 mg	8E-S012-EGB	2 Plates/Box	

#### Ordering Information

Format	Sorbent Mass	Part Number	Unit	Price
<b>Tube</b>				
	100 mg	8B-S009-EAK	1 mL (100/box)	
	200 mg	8B-S009-FBJ	3 mL (50/box)	
	500 mg	8B-S009-HBJ**	3 mL (50/box)	
	500 mg	8B-S009-HCH	6 mL (30/box)	
	1 g	8B-S009-JCH	6 mL (30/box)	

#### Giga Tube

	500 mg	8B-S009-HDG	12 mL (20/box)	
	2 g	8B-S009-KDG	12 mL (20/box)	
	5 g	8B-S009-LEG	20 mL (20/box)	
	10 g	8B-S009-MFF	60 mL (16/box)	
	20 g	8B-S009-VFF	60 mL (16/box)	


#### 96-Well Plate

	25 mg	8E-S009-CGB	2 Plates/Box	
	50 mg	8E-S009-DGB	2 Plates/Box	
	100 mg	8E-S009-EGB	2 Plates/Box	

\*\*Tab-less tubes available. Contact Phenomenex for details.

### Cyano (CN)

A polar phase with slight hydrophobic selectivity in reversed phase mode and moderate polar selectivity in normal phase mode.



#### Ordering Information

Format	Sorbent Mass	Part Number	Unit	Price
<b>Tube</b>				
	100 mg	8B-S007-EAK	1 mL (100/box)	
	200 mg	8B-S007-FBJ	3 mL (50/box)	
	500 mg	8B-S007-HBJ	3 mL (50/box)	
	500 mg	8B-S007-HCH	6 mL (30/box)	
	1 g	8B-S007-JCH	6 mL (30/box)	

#### Giga Tube

	2 g	8B-S007-KDG	12 mL (20/box)	
--	-----	-------------	----------------	--

#### 96-Well Plate

	50 mg	8E-S007-DGB	2 Plates/Box	
	100 mg	8E-S007-EGB	2 Plates/Box	