

Develosil®

By Nomura Chemical Co.

- **RP-Aqueous and RP-Aqueous-AR (C30) work with 100% aqueous mobile phases and give added polar compound retention**
- **Cyano-UG is extremely pH stable**
- **Combi-RP (C30) columns are specially packed for high flow rate Combinatorial Chemistry Applications**
- **C18, C8, Phenyl and Cyano-UG Series provide excellent peak shape for acidic, basic and chelating compounds**

Supplied worldwide (except Japan) by Phenomenex, Develosil® HPLC columns provide high efficiency, high resolution separations.

Several different C18 chemistries enable chromatographers to choose the specific column for their application requirements. Each phase provides a slightly different selectivity to help fine-tune difficult separations.

Ordering Information

Develosil Columns

Part No.	Mfr. No.	Description	Size (mm)	Price
3 µm ODS-UG (C18)				
CHO-3506	UG11346050W	ODS-UG-3	50 x 4.6	
CHO-3407	UG11346100W	ODS-UG-3	100 x 4.6	
CHO-3408	UG11346150W	ODS-UG-3	150 x 4.6	
5 µm ODS-UG (C18)				
CHO-3379	UG11546150W	ODS-UG-5	150 x 4.6	
CHO-3380	UG11546250W	ODS-UG-5	250 x 4.6	
CHO-3384	UG11540010W	ODS-UG-5 (Guard Column)	10 x 4.0	
5 µm C8-UG (C8)				
CHO-4221	UG12546150W	C8-UG-5	150 x 4.6	
CHO-4222	UG12546250W	C8-UG-5	250 x 4.6	
5 µm CN-UG (Cyano)				
CHO-4245	UG16546150W	CN-UG-5	150 x 4.6	
CHO-4246	UG16546250W	CN-UG-5	250 x 4.6	
CHO-4250	UG16540010W	CN-UG-5 (Guard Column)	10 x 4.0	
5 µm TMS-UG (C1)				
CHO-4228	UG14546100W	TMS-UG-5	100 x 4.6	
CHO-4229	UG14546150W	TMS-UG-5	150 x 4.6	
CHO-4230	UG14546250W	TMS-UG-5	250 x 4.6	
CHO-4234	UG14540010W	TMS-UG-5 (Guard Column)	10 x 4.0	
5 µm ODS-MG (C18)				
CHO-4469	MG11546150W	ODS-MG-5	150 x 4.6	
CHO-4470	MG11546250W	ODS-MG-5	250 x 4.6	
CHO-4477	MG11540010W	ODS-MG-5 (Guard Column)	10 x 4.0	
3 µm ODS-HG (C18)				
CHO-3411	HG11346050W	ODS-HG-3	50 x 4.6	
CHO-3413	HG11346100W	ODS-HG-3	100 x 4.6	
CHO-3414	HG11346150W	ODS-HG-3	150 x 4.6	
5 µm ODS-HG (C18)				
CHO-3387	HG11546150W	ODS-HG-5	150 x 4.6	
CHO-3388	HG11546250W	ODS-HG-5	250 x 4.6	
3 µm COMBI-RP				
CHO-5896	COMB346050W	Combi-RP-3	50 x 4.6	
5 µm COMBI-RP				
CHO-5902	COMB546050W	Combi-RP-5	50 x 4.6	
CHO-5904	COMB5P2050W	Combi-RP-5	50 x 20.0	

Develosil Columns (continued)

Part No.	Mfr. No.	Description	Size (mm)	Price
3 µm RP-Aqueous (C30)				
CHO-5991	RPAQ310100W	RP-AQUEOUS-3	100 x 1.0	
CHO-5992	RPAQ310150W	RP-AQUEOUS-3	150 x 1.0	
CHO-5996	RPAQ320050W	RP-AQUEOUS-3	50 x 2.0	
CHO-5998	RPAQ320150W	RP-AQUEOUS-3	150 x 2.0	
CHO-5999	RPAQ330050W	RP-AQUEOUS-3	50 x 3.0	
CHO-6001	RPAQ330150W	RP-AQUEOUS-3	150 x 3.0	
CHO-6003	RPAQ346050W	RP-AQUEOUS-3	50 x 4.6	
CHO-6006	RPAQ346150W	RP-AQUEOUS-3	150 x 4.6	
5 µm RP-Aqueous (C30)				
CHO-5683	RPAQ520050W	RP-AQUEOUS-5	50 x 2.0	
CHO-5685	RPAQ520150W	RP-AQUEOUS-5	150 x 2.0	
CHO-5686	RPAQ520250W	RP-AQUEOUS-5	250 x 2.0	
CHO-5687	RPAQ546050W	RP-AQUEOUS-5	50 x 4.6	
CHO-5689	RPAQ546150W	RP-AQUEOUS-5	150 x 4.6	
CHO-5690	RPAQ546250W	RP-AQUEOUS-5	250 x 4.6	
CHO-5695	RPAQ5P1250W	RP-AQUEOUS-5	250 x 10.0	
CHO-5699	RPAQ5P2050W	RP-AQUEOUS-5	50 x 20.0	
CHO-5696	RPAQ5P2250W	RP-AQUEOUS-5	250 x 20.0	
CHO-5672	UG17540010W	RP-AQUEOUS (Guard Column)	10 x 4.0	
3 µm RP-Aqueous-AR (C30)				
CHO-7393	RPAR310050W	RP-AQUEOUS-AR-3	50 x 1.0	
CHO-7395	RPAR320150W	RP-AQUEOUS-AR-3	150 x 2.0	
CHO-7409	RPAR346050W	RP-AQUEOUS-AR-3	50 x 4.6	
CHO-7396	RPAR346100W	RP-AQUEOUS-AR-3	100 x 4.6	
CHO-7397	RPAR346150W	RP-AQUEOUS-AR-3	150 x 4.6	
5 µm RP-Aqueous-AR (C30)				
CHO-7401	RPAR520150W	RP-AQUEOUS-AR-5	150 x 2.0	
CHO-7403	RPAR546150W	RP-AQUEOUS-AR-5	150 x 4.6	
CHO-7404	RPAR546250W	RP-AQUEOUS-AR-5	250 x 4.6	
CHO-7405	RPAR540010W	RP-AQUEOUS-AR-5	10 x 4.0	
5 µm 300 ODS-HG (Widopore C18)				
CHO-3396	HG31546250W	300ODS-HG-5	250 x 4.6	
Guard Cartridges for UG and HG columns				
CHO-3417	UG11540010C	ODS-UG-5 (4/pk)	10 x 4.0	
Guard Cartridge Holder for UG, HG, and RP-Aqueous columns				
CHO-3419	H000040010C	Guard Cartridge Holder	10 x 4.0	