

Таблица соответствия сорбентов для ТФЭ по фирмам производителям

Phenomenex <i>Strata</i>	Waters <i>Set-Pak</i>	Varian <i>Bond Elut</i>	Supelco <i>Discovery</i>	UCT	JT BAKER <i>Bakerbond</i>	IST <i>Isolute</i>	Macherey- Nagel <i>Chromabond</i>
обращенно-фазовый механизм удерживания							
C18-E	tC18	C18	DSC-18	C18	Octadecyl	C18(EC)	C18 ec
C18-U		C18 OH			PolarPlus Octadecyl Lightload	C18	C18 C18 PAH
C18-T		C18 EWP	DSC-18Lt				
C8	C8	C8	DSC-8	C8	Octyl	C8(EC)	C8
Phenyl (PH)		PH	DSC-PH	Phenyl	Phenyl	PH(EC)	
SDB-L	Porapak RDX	ENV LMS	DSC- PS/DVB		H ₂ O-phobic DVB	101	HR-P PS-RP
Strata-X	Oasis HLB	NEXUS PPL FOCUS PLEXA		StyreScreen DVB	H ₂ O-philic DVB	ENV+ ABN	Easy HR-X
нормально-фазовый механизм удерживания							
Si-1(Silica)	Silica	Silica	DSC-Silica	Silica	Silica Gel	Silica	SiOH
FL-RP	Florisil	Florisil	Florisil	Florisil	Florisil	FL	Florisil
NH₂	NH ₂	NH ₂	DSC-NH ₂	Amino Propyl	Amino	NH ₂	NH ₂
CN	CN	CN-U	DSC-CN	CN	Cyano	CN	CN
механизм удерживания по ионному обмену							
SAX		SAX	DSC-SAX	Quarternary Amine	Quarternary Amine	SAX	SB
SCX		SCX	DSC-SCX	Benzene Sulfonic Acid	Aromatic Sulfonic Acid	SCX	SA
Strata-X-C	Oasis MCX			StyreScreen DBX		CX	
Strata-X-CW	Oasis WCX						
Strata-X-AW	Oasis WAX						
WCX		CBA	DSC-WCX	Carboxylic Acid	Carboxylic Acid	CBA	PCA
Screen-C		Cetrify I		Clean Screen DAU	Narc-2	HGX	Drug
Screen-C GF		Cetrify HF		Xtract DAU			
Screen-A		Cerify II		Clean Screen DAU	Narc-1	HGX	